


PERFORMANCE. PROTECTION. EFFICIENCY.

DAR™ Filters



Filters offer protection by removing bacteria and viruses before they enter a patient's airway. Without them, a patient can be at greater risk for developing a healthcare-associated infection.

Equally important, filters reduce the number of pathogens a patient exhales into the air. This helps protect everyone, including staff, visitors, and other patients — as well as equipment.

Ventilator filters can be either electrostatic or mechanical:

- Electrostatic filters use positive and negative charges to attract and capture particles.
- Mechanical filters feature a multilayered, pleated medium that provides greater efficiency than electrostatic filters.¹

The large DAR™ mechanical filters' pleated medium significantly increases bacterial filtration efficiency,² reaching an NaCl efficiency of greater than 99.97%.³

ELECTROSTATIC FILTERS



Electrostatic filter, large

MECHANICAL FILTERS



Mechanical filter, large



Electrostatic filter, small



Mechanical filter, compact



Electrostatic filter, small, angled port



Mechanical filter, small

ELECTROSTATIC FILTERS

| | Large | Small | Small, angled port |
|---|--|-------------------------|-------------------------|
| Catalog number | 350U5865 (Without end-tidal CO ₂ sampling port) | 350U5879 | 350U19006 |
| Quantity/box | 50 | 50 | 50 |
| Recommended tidal volume | 300–1500 mL | 150–1200 mL | 150–1200 mL |
| Resistance to flow before use (ISO 9360) | | | |
| 30 L/min | 0.6 cm H ₂ O | 0.8 cm H ₂ O | 0.9 cm H ₂ O |
| 60 L/min | 1.5 cm H ₂ O | 2.1 cm H ₂ O | 2.3 cm H ₂ O |
| 90 L/min | 2.6 cm H ₂ O | 3.7 cm H ₂ O | 4.3 cm H ₂ O |
| Filtration efficiency | | | |
| Bacterial | ≥99.9999% | ≥99.9999% | ≥99.9999% |
| Viral | ≥99.999% | ≥99.999% | ≥99.999% |
| NaCl | ≥99.592%* | ≥98.096% | ≥98.096% |
| Internal volume | 99 mL | 36 mL | 44 mL |
| Weight | 35 g | 19 g | 21 g |
| Type of filtration | Electrostatic | Electrostatic | Electrostatic |

MECHANICAL FILTERS

| | Small | Compact | Large | Large w/o gas sampling port |
|---|-------------------------|-------------------------|-------------------------|--------------------------------|
| Catalog number | 351U5979 | 351U5878 | 351U5410 | 351U5856 |
| Quantity/box | 50 | 50 | 50 | 50 |
| Recommended tidal volume | 150–1200 mL | 200–1500 mL | 300–1500 mL | 300–1500 mL |
| Resistance to flow before use (ISO 9360) | | | | |
| 30 L/min | 1.2 cm H ₂ O | 0.8 cm H ₂ O | 0.8 cm H ₂ O | 0.8 cm H ₂ O |
| 60 L/min | 2.7 cm H ₂ O | 1.9 cm H ₂ O | 2.0 cm H ₂ O | 2.0 cm H ₂ O |
| 90 L/min | 4.5 cm H ₂ O | 3.2 cm H ₂ O | 3.6 cm H ₂ O | 3.6 cm H ₂ O |
| Filtration efficiency | | | | |
| Bacterial | ≥99.9999% | ≥99.9999% | ≥99.9999% | ≥99.9999% |
| Viral | ≥99.999% | ≥99.9999% | ≥99.9999% | ≥99.9999% |
| NaCl | ≥99.512%* | ≥99.747% ³ | ≥99.978%* | ≥99.978%* |
| Internal volume | 42 mL | 66 mL | 92 mL | 92 mL |
| Weight (approx.) | 24 g | 39 g | 47 g | 47 g |
| Type of filtration | Mechanical | Mechanical | Mechanical | Mechanical |

*Internal testing Mirandola (various 2005–2008).

1. Cann C, Hampson MA, Wilkes AR, Hall JE. The pressure required to force liquid through breathing system filters. *Anaesthesia*. 2006;61(5):492–497.
2. Wilkes AR. Measuring the filtration performance of breathing system filters using sodium chloride particles. *Anaesthesia*. 2002;57(2):162–168.
3. Nelson Laboratories Inc. Sodium chloride aerosol testing of breathing system filters (BSF). Lab No. 399951A. 1 Amended. January 2008.

© 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 09/2016-16-AW-0030-[WF#1065275]