

# Dale® BreezeLock® Endotracheal Tube Holder



For Comfortable & Secure Stabilization

# Dale® Medical Products



**Dale Medical Products, Inc.**  
40 Kenwood Circle • Suite 7  
Franklin, MA 02038-3298, USA  
T: +1-508-695-9316  
[www.dalemed.com](http://www.dalemed.com)

EC REP

**MDSS GmbH**  
Schiffgraben 41  
D-30175 Hannover, Germany



CH REP

**MDSS CH GmbH**  
Laurenzenvorstadt 61  
5000 Aarau, Switzerland

(Croatian, Danish, Dutch, English, Finnish, French, German, Greek, Italian, Japanese, Korean, Norwegian, Polish, Portuguese, Spanish, Swedish, Thai, Turkish)

# Dale® Medical Products

## BreezeLock® Endotracheal Tube Holder

fits endotracheal tube sizes 7.0-10.0mm

**Contents:** 1 of each - neckband; 273 Adhesive Base; roll of tape

### Application Instructions

- Clean tube. Let dry.
- Apply Adhesive Base. **(A)** Follow directions on 273 package.
- Press neckband onto hook portion of the Adhesive Base, making sure the TRACK is centered above the patients mouth. **(B)**
- Open the fastener tabs located at both ends of the neckband and wrap the band around the head just below the ears. Secure the fastener tabs to the outside of the neckband.
- Lift DALE OVERLAY. Position the CHANNEL over the ET Tube. Press DALE OVERLAY across TRACK and Neckband to secure into position.
- Steps to tape the tube to the CHANNEL:
  - With the tube resting freely below the CHANNEL, wrap 6" long piece of provided adhesive tape 1x around the tube. **(C)**
  - Press TUBING CHANNEL down onto the tube.
  - Firmly wrap tape in a forward spiral motion around tube and CHANNEL. **(D)**
  - Optional: Wrap tape a minimum of 1" beyond the end of the TUBING CHANNEL.
- Ensure neckband security by readjusting the fastener tabs as necessary.
  - Change tape on CHANNEL as necessary.
  - Change device after 3 days or sooner if necessary.
  - To prevent skin ulcerations or skin tears, assess patient skin and medical conditions prior to applying. Monitor daily or more frequently.

### Tube Repositioning

While holding the taped CHANNEL with one hand, lift DALE OVERLAY off TRACK. Slide CHANNEL to move tube to desired location. Press DALE OVERLAY back down onto TRACK and neckband to secure the tube's position in place.

### Removal

Device should be removed with care to avoid skin trauma. Re-apply skin prep if another device is to be applied.

Dale **recommends repositioning the endotracheal tube often** to help prevent injury to the lips and underlying tissues due to unrelieved pressure.

Dale **strongly recommends supporting the ventilator tubing** to reduce pressure on the patients face.

For single patient use. Disposable.

For an application video please visit [www.dalemed.com](http://www.dalemed.com)

Not made with natural rubber latex.

Patent [www.dalemed.com/patents](http://www.dalemed.com/patents)

Dale & BreezeLock are registered trademarks of Dale Medical Products, Inc. in the USA, EU, and UK.

