

Cutimed® DebriClean

Fast and effective wound debridement


**> 99% biofilm
removal with
just 4 wipes^{1,2}**



Meet Cutimed® DebriClean

Fast and effective wound debridement

In chronic wounds, devitalised tissue and biofilm forms a physical barrier that inhibits the formation of new tissue. For the wound to heal properly, the devitalised tissue and biofilm must be removed.³ That's why every layer of the **Cutimed® DebriClean** pad has been specially designed to support you in the mechanical debridement process and lead to **more effective wound healing**.




* Product is compatible with commonly used wound cleansing solutions based on PHMB, sodium hypochlorite or hypochlorous acids, octenidine, 0.9% medical saline solution, water and Ringer's solution



Cutimed® DeбриClean


all-purpose debridement

Thorough wound cleansing is the foundation for successful wound healing. **Cutimed® DeбриClean** sets the stage for the wound healing process with its unique and innovative design. Due to its combination of gentle looped monofilament fibres and abrasive fibres, bacteria is absorbed and firm slough removed, making **Cutimed® DeбриClean** a suitable wound debridement pad for wound bed preparation.




Effective bacterial binding

In *in vitro* testing, **Cutimed® DebrisClean** excelled in bacterial binding when compared to standard debridement options.




How many bacteria remained bound to the debridement pad after rinsing?



*Initial concentration 2.5 Mio cfu/ml *Pseudomonas aeruginosa* DSM1117. Contact time 15 min. intensive rinsing with Tween 20 (non ionic tenside)

Strong cleaning efficiency

When tested *in vitro* for efficiency, **Cutimed® DebrisClean** provided very good results compared to standard debridement options. The cleaning efficiency of each product was tested through the mechanical wiping of protein crust, protein networks and hyperkeratosis.



Choose Cutimed® DebricClean for fast and effective wound debridement



Strong cleaning efficiency^{1,2}



Effective bacterial binding⁴



> 99% biofilm removal²



Ordering information

Code	Size	Pack Qty	Shipping Qty	NHS code
76170-00	9 x 13 cm	1 x 5	12 boxes	415-3631

¹ JB4 40031810, In vitro evaluation of the cleansing effect of Cutimed® DebricClean and Cutisoft Cotton. August 2018, BSN medical data on file.

² JB4 40031908, In vitro evaluation of removal of viable biofilm by Cutimed® DebricClean and Cutisoft Cotton. August 2018, BSN medical data on file.

³ EWMA, EWMA document: Debridement. 2013, European Wound Management Association. p. 1-49.

⁴ JB4 40031811, In vitro evaluation of the bacterial adhesion to Cutimed® DebricClean and Cutisoft® Cotton. July 2018, BSN medical data on file.