

NOL™ (Nociception Level) Technology

Intraoperative Pain Monitoring
for Optimizing Analgesia



Accurate & Objective
assessment of
nociception



Optimization & Personalization
of analgesic treatment



Faster Response Time
through earlier detection
of nociception



Standardization
among physicians

THE CHALLENGE

Assessing Intraoperative Pain Reliably

Treating pain is at the heart of medicine. It is an essential role of every clinician, since unmanaged pain can delay recovery, increase morbidity and mortality, and overburden healthcare resources.

Nociception - Pain in Anaesthetized Patients


During general anaesthesia, a patient's body reacts to painful stimuli - although it is not consciously recognized. This intraoperative pain can stress the body¹⁻⁴ and worsen pain after surgery^{4,5}. As the patient can't communicate - it is hard for clinicians to evaluate.

Missing a Piece?

While Hypnosis and Immobility are continuously and specifically monitored, analgesia is assessed indirectly through changes in hemodynamic and clinical parameters (Heart Rate, Blood Pressure, Sweating, Tearing, etc.).

These parameters are not specific to nociception, and may change in response to other causes as well.

Consequently, patients may be given insufficient analgesia, which may promote postoperative pain^{4,5}, or excessive analgesia, which may cause respiratory complications^{1,6}, nausea and vomiting⁷, hyperalgesia^{8,9}, and other complications¹⁻¹¹.



” Getting the right dose of anti-nociceptive medications matters. Too little, and patients wake up in pain. Too much, and patients are at risk of drug-related complications. ”

Dr. Daniel Sessler, Head of the Department of Outcomes Research, Cleveland Clinic, Ohio, USA. A member of Medasense's advisory board.

Each Year, Worldwide:



UP TO **50%**

of surgical patients suffer from moderate to severe post-operative pain¹²⁻¹³.



12%

of surgical patients suffer from adverse events due to analgesic medications, leading to¹³:

3.3 extra days of hospitalization

27% extra cost per patient

Increase in re-admissions



Monitoring the intraoperative nociception level may enable clinicians to personalize and optimize the analgesic administration, thus avoiding excessive use or under use of opioids, that may result in significant complications.

THE SOLUTION

NOL™ (Nociception Level) Technology

Multi-Parametric Approach

Recognizing the complex nature of pain, the NOL™ technology considers multiple nociception-related physiological parameters and its various derivatives.



Heart Rate



Heart Rate Variability



Skin Conductance Level



Skin Conductance Fluctuations



Pulse Wave Amplitude



Peripheral Temperature



Motion



Parameters Derivatives

Advanced Algorithms

Advanced algorithms process the multiple data streams, identify the pain-related patterns, and reflect a patient's nociceptive state¹⁴⁻¹⁶.


” Multi-variable approaches appear to be superior predictors of pain intensity and intra-operative nociception to any individual parameter alone. ”

Cowen R et al., (2015) Assessing pain objectively: the use of physiological markers. *Anaesthesia* 70(7):828-47.

NOL™ (Nociception Level) Index

- An **objective indicator** for the presence and severity of nociception.
- The only **multi-parameter nociception index**, based on the physiological integrated response to noxious stimuli.
- **Calibrated to an individual's baseline**.
- **Clinically validated** as superior to other nociception indicators¹⁴⁻¹⁶.

The NOL™ index values are represented on a scale of 0 to 100



0 = Absence of Nociception → 100 = Extreme Nociception

PMD-200™

Intraoperative Nociception Monitoring

Leveraging patented NOL™ technology, the PMD-200™ is a non-invasive and continuous pain monitoring device.

By using the PMD-200™ in operating rooms, where patients under general anaesthesia are unable to communicate their pain, clinicians are able to assess nociception, and better titrate analgesic medications - avoiding excessive use or underuse, that may result in significant complications.



A Non-Invasive Finger Probe

The finger probe continuously acquires physiological signals through 4 sensors (Photoplethysmograph, Galvanic Skin Response, Temperature, and Accelerometer) and transmits them to the PMD-200™ for advanced processing.

PMD-200™ Bedside Monitor

Pattern recognition algorithms analyze dozens of measurable changes in the nociception-related physiological parameters, and reflect the information in a single nociception level index - the NOL™.

” I truly believe in the NOL™ technology. I am convinced it will improve the quality of our anaesthesia and impact on patients’ recovery.”

Dr. Philippe Richebé, Professor and Anaesthesiologist at Maisonneuve-Rosemont Hospital, University of Montreal, Quebec, Canada. Principal Investigator. A member of Medasense’s advisory board.

CLINICAL EVIDENCE

Superior Nociception Assessment

Multiple studies conducted worldwide and key results published in peer reviewed medical journals, have validated the NOL™ Index to:

- Discriminate, with high sensitivity and specificity, between noxious and non-noxious stimuli¹⁴⁻¹⁶.
- Grade different levels of nociception¹⁴⁻¹⁶.
- Correlate with the analgesic state of the patient^{15,17}.
- Demonstrate superiority over commonly used nociception indicators¹⁴⁻¹⁶.

“ Our study confirmed that the NOL™ gives a better assessment of the patient’s nociceptive state than the parameters we currently use in the operating theatre. ”

Dr. Ruth Edry, Department of Anaesthesiology, Rambam Medical Centre, Haifa, Israel.





Medasense’s updated publications, abstracts and presentations can be found at: www.medasense.com/clinical-data/



Superiority of the NOL™ index to detect and discriminate between various noxious stimuli, compared to commonly used parameters¹⁵

(N=71; ASA I - III; Ages 18-80; BIS target 45+/-5 ; Elective surgery under general anaesthesia)

Only NOL™ correctly scored the level of nociceptive reaction with:
non-noxious stimulus NOL™ < incision NOL™ < intubation NOL™ (p < 0.05).




○ Prior to stimulation ● Following stimulation # Unpaired * Paired

The NOL™ index reliably measures the changes in the nociceptive response at different remifentanyl concentrations (inter-patients)¹⁵

(N=71; ASA I - III; Ages 18-80; BIS target 45+/-5 ; Elective surgery under general anaesthesia)

The NOL™ index remains unaffected under non-noxious conditions, regardless of remifentanyl concentration and decreases for the same noxious stimulus with increasing remifentanyl concentrations.


This analysis indicates that the NOL™ is a reliable measure of nociception and is not affected by the hemodynamic effects of remifentanyl.



The NOL™ index correlates with increased dosage of analgesics (intra-patient)¹⁷

(N=40; ASA I - III; Age >18; Elective abdominal surgery under general anaesthesia and epidural analgesia)

The magnitude of the NOL™ index response to standardized nociceptive stimulus, decreases with higher doses of remifentanyl.




The NOL™ index outperforms commonly used parameters (HR, PPGA) and SPI (GE)¹⁴

(N=58; ASA I - III; Ages 18-75; Entropy target <60 ; Elective surgery under general anaesthesia)

In this ROC analysis NOL™ outperforms other parameters and indices to discriminate between noxious and non-noxious stimuli.

AUC for NOL™ absolute values were the highest: 0.93.



NOL™ Technology

Setting a New Standard in Nociception Assessment

” For the first time we are able to titrate analgesic medications to patients’ needs. ”

” It is clear that current pain assessment approaches in the operating room are limited. The NOL™ index provides anaesthesiologists with a decision support tool to objectively assess and optimize the treatment of a patient’s nociception pain. ”

Prof. Albert Dahan, MD, PhD, Department of Anaesthesiology, Leiden University Medical Center, The Netherlands. Principal Investigator.

References

1. Desborough JP, (2000), The stress response to trauma and surgery, Br J Anaesth. 85:109-117.
2. Dickerson DM, (2014), Acute pain management. Anesthesiol Clin. 32:495-504.
3. McGreevy K, (2011), Preventing Chronic Pain following Acute Pain: Risk Factors, Preventive Strategies, and their Efficacy. Eur J Pain Suppl. 5:365-337.
4. Kehlet H et al., (2006), Persistent postsurgical pain: risk factors and prevention, Lancet 367(9522):1618-25.
5. Kehlet H et al., (1997), Multimodal approach to control postoperative pathophysiology and rehabilitation. Br J Anaesth 78:606-617.
6. Lee LA et al., (2015), Postoperative Opioid-induced Respiratory Depression: A Closed Claims Analysis. Anesthesiology 122:659-665.
7. Smith HS et al., (2014), Opioid induced nausea and vomiting. Eur J Pharmacol. 722:67-78.
8. Joly V et al., (2005), Remifentanil-induced postoperative hyperalgesia and its prevention with small dose ketamine. Anesthesiology 103:147-55.
9. Koppert W et al., (2007), The impact of opioid-induced hyperalgesia for postoperative pain. Best Practice & Research Clinical Anaesthesiology Volume 21, Issue 1, Pages 65-83.
10. Gan TJ et al., (2014), Incidence, patient satisfaction, and perceptions of post-surgical pain: Results from a US national survey. Current Medical Research and Opinion. 30:149-6.
11. Benyamin R et al., (2008), Opioid complications and side effects. Pain Physician. 11:S105-120.
12. Pogatzki-Zahn E et al., (2015), A Prospective Multicenter Study to Improve Postoperative Pain: Identification of Potentialities and Problems. PLoS ONE 10 (11): e0143508.
13. Oderda G et al., (2013), Effect of Opioid-Related Adverse Events on Outcomes in Selected Surgical Patients. J Pain Palliat Care Pharmacother, 27:62-70.
14. Edry R et al., (2016), Preliminary Intraoperative Validation of the Nociception Level Index: A Noninvasive Nociception Monitor. Anesthesiology 125:193-203.
15. Martini C et al., (2015), Ability of the Nociception Level (NOL), a multiparameter composite of autonomic signals, to detect noxious stimuli during propofol-remifentanil anesthesia. Anesthesiology. 123:524-534.
16. Ben-Israel N et al., (2013), Monitoring the nociception level: a multi-parameter approach. J Clin Monit Comput. 27(6):659-68.
17. Richebé P et al., (2016), Nociception Level (NOL) index alteration after standardized nociceptive stimulus decreases with higher doses of remifentanil. ESA Abstract (European Society of Anaesthesiology), London, UK, 2016.



info@medasense.com | www.medasense.com

Medasense Biometrics Ltd. develops innovative systems and applications to objectively assess the physiological response to pain. Objective assessment of pain may personalize pain management, leading to optimized care and improved clinical outcomes.