sinus-SuperFlex-DS · sinus-Repo-DS

Self-expanding nitinol stent system

Solutions for the
Hypoplastic Left
Heart Syndrome (HLHS)

DUCTAL STENTING

sinus- DS


The Hypoplastic Left Heart Syndrome

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The Hypoplastic Left Heart Syndrome

Hypoplastic Left Heart Syndrome (HLHS) is the term used for a variety of simultaneously occurring malformations of the heart and aorta in neonates. The left ventricle of the heart is underdeveloped and completely closed off from the aorta, so that blood flow to the aorta is provided by the ductus arteriosus to ensure that essential organs are supplied with blood. HLHS is among the most serious congenital heart defects and is fatal in 99% of cases if left untreated.

A number of options are available for successful treatment of HLHS. The hybrid procedure reduces the patient stress to a minimum. A heart-lung machine is not required and the number, complexity and duration of interventions can be reduced.

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Exceptional solutions for exceptional situations

sinus-DS stents are designed specifically for treatment of Hypoplastic Left Heart Syndrome and are approved class III devices (93/42/EEC). Both versions are based on proven sinus-technology and produced seamlessly from a single piece of nitinol without any joints or welds. This process enhances resistance to stent breakage and maximises radial force of even extremely short stents. A special electro-polishing process ensures ultra-smooth surfaces and additionally reduces thrombogenesis.

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The self-expanding sinus-SuperFlex-DS stent in open-cell design is designed specifically to effect temporary patency of the ductus arteriosus. With flexibility and radial force, the stent adjusts to the curves of the vessel wall.

- Flexibility
- Radiopaque markers
- Vessel adaptation
- Electro-polishing

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The self-expanding sinus-Repo-DS stent in closed-cell design is for additional use, in case of a coarctation of the descending aortic arch above the ductus connection.

- Radial force
- Radiopaque markers
- Repositionable (up to 80% of nominal stent length)
- Electro-polishing


sinus-SuperFlex-DS

Originalgröße Original size

sinus-Repo-DS

Originalgröße Original size

The first ductal stent system with class III certification worldwide (93/42/EEC)

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The procedure

The hybrid procedure is a cooperative procedure combining techniques of catheter intervention and cardiac surgery.

After birth, both branches of the pulmonary artery are surgically narrowed using bands, reducing the high flow of blood to the lungs and directing the blood through the patent ductus arteriosus into the aorta. A stent is implanted to maintain blood flow through the ductus arteriosus.

The open-cell design of the sinus-SuperFlex-DS stent provides a very good combination of flexibility and radial force to establish temporary patency in the ductus arteriosus.

In cases of coarctation of the descending aortic arch where stenting of the ductus arteriosus is required, the sinus-Repo-DS stent is used in addition. The closed-cell design of this stent maximises radial force and enables the repositioning of the stent (up to 80% of the nominal length).

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Maximum flexibility

The highly flexible 4F application system of the sinus-DS stents provides three implantation options:

transvenous : over the femoral vein

transarterial : over the femoral artery

transpulmonal : during hybrid surgical procedure, after application of bilateral banding of the pulmonary artery

The application system features an exceptionally easy-to-use pull-back-system. Radiopaque markers, braided sheath and atraumatic soft tip facilitate guidance and positioning. To ensure accurate positioning, the stent is fixed in the application device by the anti-jump system - a integrated mechanism which prevents the stent from jumping.

- Pull-back-system
- Braided outer sheath
- Anti-jump-technique
- Atraumatic soft tip

sinus-SuperFlex-DS

For temporary patency of the ductus arteriosus

Bestell-Nr. order code référence	∅ Stent	↔ Stent	System
	mm	mm	F / cm
8804-2015	4	15	4F / 85 cm
8804-2018	4	18	4F / 85 cm
8804-2020	4	20	4F / 85 cm
8804-2022	4	22	4F / 85 cm
8804-2024	4	24	4F / 85 cm
8805-2015	5	15	4F / 85 cm
8805-2020	5	20	4F / 85 cm
8805-2024	5	24	4F / 85 cm
8806-2015	6	15	4F / 85 cm
8806-2020	6	20	4F / 85 cm
8806-2024	6	24	4F / 85 cm
8807-2012	7	12	4F / 85 cm
8807-2015	7	15	4F / 85 cm
8807-2018	7	18	4F / 85 cm
8807-2020	7	20	4F / 85 cm
8807-2022	7	22	4F / 85 cm
8807-2024	7	24	4F / 85 cm
8808-2012	8	12	4F / 85 cm
8808-2015	8	15	4F / 85 cm
8808-2018	8	18	4F / 85 cm
8808-2020	8	20	4F / 85 cm
8808-2022	8	22	4F / 85 cm
8808-2024	8	24	4F / 85 cm
8809-2015	9	15	4F / 85 cm
8809-2020	9	20	4F / 85 cm
8809-2022	9	22	4F / 85 cm
8809-2024	9	24	4F / 85 cm

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Verpackungseinheit: 1 Stück
conditionnement: 1 unité
box: 1 unit

sinus-Repo-DS

For treatment of existing coarctations

Bestell-Nr. order code référence	∅ Stent	↔ Stent	System
	mm	mm	F / cm
8905-2009	5	9	4F / 85 cm
8906-2009	6	9	4F / 85 cm

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