

Worldwide Addresses

Germany

Krablerstrasse 127
D - 45326 Essen
Tel: (+49) 201 861 860 - 0
Fax: (+49) 201 861 86 - 12
office@de.gbo.com

Switzerland

St. Leonhardstrasse 31
9000 St. Gallen
Tel: (+41) 71 228 55 22
Fax: (+41) 71 228 55 21
office@ch.gbo.com

United Kingdom

Stroudwater Business Park
Brunel Way
Stonehouse, Glos. GL10 3SX
Tel: (+44) 14 53 82 52 55
Fax: (+44) 1453 82 72 77
office@uk.gbo.com

Spain

Parque Empresarial INBISA-1
Av. Somesierra, 22-nave 1A
San Sebastian de los Reyes
28700 Madrid
Tel: (+34) 91 652 77 07
Fax: (+34) 91 652 33 35
info@vacuette.es

The Netherlands

Albert Einsteinweg 16
2408 AR Alphen a/d Rijn
Tel: (+31) 172 42 09 00
Fax: (+31) 172 44 38 01
office@nl.gbo.com

Brazil

Rua Affonso Pansan, 1967 -
13.473 - 620 Vila Bertini
Americana - Sao Paulo
Tel : (+55) 19 3468 9600
Fax : (+55) 19 3468 9601
info@br.gbo.com

U.S.A

4238 Capital Drive
Monroe, North Carolina 28112
Tel: (+1) 704 261 7800
Fax: (+1) 704 261 7899
Customer service: 888-286-3883
office@us.gbo.com

Hungary

Fertősor 7
H – 9200 Mosonmagyaróvár
Tel: (+36) 96 213 088
Fax: (+36) 96 213 198
office@hu.gbo.com

France

3 - 7 avenue du Cap Horn
BP 31 - Les Ulis
91941 Courtaboeuf
Tel: (+33) 1 69 86 25 25
Fax: (+33) 1 69 86 25 35
office@fr.gbo.com

Portugal

R/C Bloco A
Areia - Arvore
Rua de Estrada Velha, 1281
Tel: (+351) 252 647 721
Fax: (+351)252 647 722
info@vacuette.pt



VACUETTE® - "Take the Original"

For decades glass was the raw material used for the production of evacuated blood collection tubes. However, in the mid 1980's, Greiner Bio-One recognised the safety risks for medical personnel that are associated with the use of glass, and was the **first producer** with the technical competence to make an innovative specimen collection system manufactured from PET, a virtually unbreakable highly transparent plastic.



VACUETTE®

Blood Collection System

Innovation for
Your Safety!

See for yourself the advantages that "The ORIGINAL" can offer:

- ⊕ Worldwide technological leadership in the production of specimen collection systems made of plastic
- ⊕ Special thick-walled tube design to guarantee a longer shelf-life
- ⊕ Innovative safety products to protect your health
- ⊕ Complete product range made of virtually unbreakable plastic
- ⊕ Flexibility for custom-made solutions with highest quality standards

"ONE STEP AHEAD"

VACUETTE® users profit from one of the most modern and highly efficient production facilities in the world. The continual development of new technologies and the state of the art production facilities demonstrate the innovative strengths of the company.

An ingenious quality management system guarantees the high standards of Greiner Bio-One. The use of VACUETTE® products ensures the safety of the user as well as the health and well-being of a magnitude of healthcare workers and patients.



Austrian headquarters of Greiner Bio One

Preparing for Blood Collection

2

a) Patient identification

This is performed by a comparison between the patient's test order form and the patient identification number, barcode, wristband number or other objective criteria.

b) Position

The patient should be suitably positioned for venipuncture (either sitting or recumbent), the position should be maintained for a minimum of 15 minutes prior to performing the venipuncture.

c) Preparation of the Collection Material

Prior to performing the venipuncture, the following items must be prepared:

- **VACUETTE®** Blood Collection System (consists of **VACUETTE®** Multiple Use Drawing Needle, Standard Tube Holder, **VACUETTE®** Blood Collection Tubes)
- Sterile disposable gloves
- Sterile swab
- Disinfectant or alcohol solution
- Adhesive bandages
- **VACUETTE®** Tourniquet
- **VACUETTE®** Sharps Disposal Container
- Label for patient identification (Timing of labelling varies from country to country)




Selection of a Puncture Site

3

► Priority list:

- 1) Blood collection from the antecubital area of the arm (Median, Basilic or Cephalic veins)



- 2) Blood collection from the dorsal side of the hand (Dorsal veins).



- 3) Blood collection from the dorsal surface of the foot (Venous arch)



Prior to making the final selection of a site for venipuncture, an inspection of the proposed area is necessary. The selection sequence should correspond to the priority list; whereby 1) and 2) are suitable in 95 % of cases and provide a satisfactory outcome.

Venipuncture

4

▶ Lightly tap the vein (only in the case of a non-prominent vein).



▶ Vein stasis with a tourniquet – maximum duration 1 minute.



▶ Disinfect the puncture site (Allow the disinfectant to thoroughly dry).



▶ **Venipuncture** – Insert the threaded needle into the holder and then insert into the patient's vein. The patient arm should be inclined in a downward position.



▶ With the second hand, the vacuum tube should be inserted into the holder (the tube cap must point upwards). Ensure that the rubber stopper is fully penetrated. Release tourniquet as soon as blood begins to flow.



Venipuncture

▶ For patients with prominent veins it is recommended to use the following **VACUETTE®** standard blood collection products:

- **VACUETTE® Multiple Use Drawing Needle** with uniquely sharpened faceting (three needle gauges available in 20, 21 and 22 G) for a patient friendly, pain free blood collection.



- **Standard Tube Holder** with economic design. The specially adjusted surface area allows for improved handling of the holder during blood collection.



- **VACUETTE® Blood Collection Tubes**, available exclusively in PET plastic. All tubes are available with a safety cap. The use of a vacuum system eliminates the possibility of back flow occurring during blood collection.



Venipuncture

► For patients with difficult veins and respectively patients presenting an increased infection risk, it is recommended to use the following **VACUETTE®** products:

- **Safety Blood Collection Set** to protect your health. The simple to use safety mechanism provides reliable protection against needlestick injuries. The Safety Blood Collection Set is available as a standard version or as two special versions (with Luer adapter and with tube holder).



- **HOLDEX® Single-Use Holder** with eccentric nozzle. This holder is especially suitable for difficult vein conditions. The puncture is guaranteed through the off-centre cannula connection providing an optimal puncture angle for a patient friendly blood collection. The **HOLDEX®** Single-Use Holder can be used together with a Luer needle or with a blood collection set.




- **VACUETTE® QUICK-SHIELD Safety Tube Holder**; the practical safety holder with integrated protective cover to prevent needlestick injuries. Activation of the protective cover takes place when the needle is removed from the patients-vein. The holder should be used exclusively with **VACUETTE®** Multiple Use Drawing Needles.



Recommended Order of Draw for Multiple Specimen Collection:

5



*Coagulation tubes may be the first tube to be drawn for routine testing only (PT and aPTT).

Note: Always follow your facility's protocol for order of draw.

Troubleshooting

6

► What should you do if no blood flows into the tube?

| Possible cause | Solution |
|--|---|
| The bevel of the needle tip is sucked against the wall of the vein | Gently rotate the needle within the lumen of the vein |
| The needle penetrated the vein wall | Gently pull both the tube holder and the needle backwards |
| The needle is not fully within the vein | Gently push the needle forwards |
| The tourniquet was too tight or in place too long | Loosen the tourniquet |
| The tube was already used, or was previously opened | Dispose of and select a new tube |

► Blood flow ceases midway through the collection

| Possible cause | Solution |
|---|--|
| The tube was removed from the holder too soon | Reinsert the tube into the holder until the vacuum is totally depleted |
| Suction is too strong for the vein (collapsed vein) | Pull the tube out of the holder for a second and then reinsert it |
| The needle position has altered during the procedure, or the needle is outside the vein | Repeat venipuncture at different site when haemotoma occurs. |

Troubleshooting


► Haemolytical sample material

| Possible cause | Solution |
|--|---|
| Too long stasis of the vein (longer than 1 minute) | Exact control of stasis time (maximum 1 minute) |
| Transfer from a syringe into a tube | For save blood transfer, use the VACUETTE® Blood Transfer Unit |
| Too intense mixing of the sample | Gently invert the tube 8 times (coagulation tubes 4 times) |
| Tubes, that are not adequately filled | Ensure that the tube is correctly filled to the fill mark on the tube label |

Mixing of Specimen Material and Labelling of Tubes

7

- Following blood collection, all tubes should be gently inverted 8 times (coagulation tubes 4 times). Thorough mixing is necessary to ensure adequate performance of the tube contents (additive) with the blood sample. A full inversion is when the air bubble moves from one end of the tube to the other.




- In order to ensure unique identification of the specimens, it is necessary to write on the tube labels or use a barcode system to label the tubes.

Transport

8

- ▶ The recommended transport and storage temperature for tubes prior to use is 4-25 °C (40-77° F). Exceeding the recommended storage temperature may lead to impairment of the tube quality.
- ▶ Avoid direct exposure to sunlight in storage and during transportation of samples, especially light-sensitive analytes such as Bilirubin.
- ▶ For safe transport, it is recommended that the **VACUETTE®** Transport Boxes - especially developed for this purpose - be used in combination with the appropriate transport carton.



Centrifugation

9


Centrifugation recommendations for **VACUETTE®** Blood Collection Tubes

| | Centrifuge speed | Time |
|----------------------------|------------------|---------|
| Coagulation tubes | | |
| Platelet rich plasma (PRP) | 150 g | 5 min. |
| Platelet poor plasma (PPP) | 1500 – 2000 g | 10 min. |
| Platelet free plasma (PFP) | 2500 – 3000 g | 20 min. |
| Serum tubes | | |
| Serum Sep C/A | min. 1500 g | 10 min. |
| Serum Beads C/A | 1800 g | 10 min. |
| Plasma tubes | | |
| Heparin Sep | 2000 - 3000 g | 15 min. |
| EDTA Sep | 2200 g | 15 min. |
| EDTA Sep | 1800 – 2200 g | 10 min. |

Serum tubes should be centrifuged within 30 minutes after blood collection.

Important:

The type of centrifuge being used can influence the composition of the gel barrier. Through the use of a swing-out rotor centrifuge in comparison to a fixed-angle centrifuge, a more solid gel barrier will be achieved. Centrifugation should be performed in a cooled centrifuge (15 – 24 °C).



Handling **VACUETTE®** Tubes with Safety Cap 10

- ▶ Opening **VACUETTE®** Blood Collection Tubes

- 1) Hold the tube firmly in one hand (use a solid base to support the arm)
- 2) Twist the safety cap with the other hand so the cap is loosened
- 3) Carefully open the tube with a gentle twist-pull motion. **VACUETTE®** PREMIUM Tubes are opened with a short twist movement.


Note: Too long storage of opened tubes can lead to evaporation and therefore false analysis results!

- ▶ Closing **VACUETTE®** Blood Collection Tubes

- 1) Place the safety cap on the tube
- 2) Press the cap onto the tube with the thumb (so it is firmly seated). **VACUETTE®** PREMIUM Tubes are closed with a twist movement.




VACUETTE® Handling of Testing Material 11




Tube Dimensions

12

75mm tubes
specimen volumes



100mm tubes
specimen volumes



VACUETTE® Safety Caps

13

► Standard cap

Snap cap



- Prevents aerosol effect
- Absolute transport security
- Simple re-capping
- Prevents contact with patient blood

Application Areas for VACUETTE® Tubes

14

| VACUETTE® tube type | Colour-coding of cap | Additive | Intendend purpose |
|---------------------|----------------------|---|---|
| Serum | | Clot Activator | Determinations in serum for clinical chemistry, microbiological serology, immunology, TDM |
| Serum Sep | | Clot Activator and Sep | Determinations in serum for clinical chemistry, microbiological serology, immunology, TDM |
| Serum Beads | | Clot Activator and Beads | Determinations in serum for clinical chemistry, microbiological serology, immunology |
| Serum Crossmatch | | Clot Activator | Determinations in serum for crossmatch testing |
| Plasma | | Sodium Heparin Lithium Heparin Ammonium Heparin | Determinations in heparinised plasma for clinical chemistry |
| Plasma Sep | | Lithium Heparin and Sep | Determinations in heparinised plasma for clinical chemistry |
| EDTA | | K ₂ EDTA K ₃ EDTA | Determinations in EDTA whole blood for haematology |
| EDTA Crossmatch | | K ₃ EDTA | Determinations in EDTA whole blood for crossmatch testing |
| EDTA Sep | | K ₃ EDTA / Sep | Determinations in EDTA plasma for molecular biological identification of viruses, parasites und bacteria |
| Coagulation | | Citrate Solution (3.2%) Citrate Solution (3.8%) | Determinations in citrated plasma for coagulation testing |
| CTAD | | CTAD (3.2%) | Determinations in citrated plasma for coagulation testing where the artificial entry of platelet factors into the plasma is avoided |
| Glucose | | Anticoagulant Glycolysis inhibitor | Determinations in stabilised anticoagulated whole blood or plasma for glucose and lactate testing |
| Trace Elements | | Clot Activator Sodium Heparin | Determinations in serum / heparinised plasma for trace elements testing |
| Blood Grouping | | ACD-A ACD-B CPDA | Determinations in ACD / CPDA whole blood for blood grouping |

VACUETTE® offers a wide variety of supportive education and training material for the sample collection procedure. Amongst others, these include:

15

- ▶ **VACUETTE®** Tube Summary Chart _____ Art. No. 980015
- ▶ **VACUETTE®** Instructions for Use Leaflet _____ Art. No. 980055
- ▶ **VACUETTE®** Hygiene Compendium _____ Art. No. 980056
- ▶ **VACUETTE®** Blood Collection Techniques Booklet _____ Art. No. 980063
- ▶ **VACUETTE®** Safety Brochure _____ Art. No. 980124
- ▶ **VACUETTE®** Preanalytics Manual _____ Art. No. 980183
- ▶ **VACUETTE®** Analyte Chart _____ Art. No. 980196
- ▶ **VACUETTE®** GBO Company Product Presentation DVD _____ Art. No. 980434



For further information, please contact your nearest **VACUETTE®** distributor.



VACUETTE®
one step ahead ▶

BLOOD COLLECTION SYSTEM

An Overview of the Essential Points

Checklist:

| | |
|---|--|
| Storage conditions in storeroom | 4 – 25°C, protected from direct sunlight |
| Expiry Date | Not exceeding end of month stated on tube and packaging |
| Blood collection | According to Instructions for Use |
| Inversion of tubes immediately after blood collection | Gently 8 times (coagulation tubes 4 times) |
| Waiting period prior to centrifugation of serum tubes | Minimum 30 min. |
| Visual control of complete coagulation of serum tubes prior to centrifugation | In certain blood samples, the clotting speed in serum tubes may be clearly delayed (i.e. anticoagulant therapy, missing coagulation factor,) the waiting period prior to centrifugation may be correspondingly delayed |
| Centrifugation | According to Instructions for Use |
| Further transport/dispatch of centrifuged gel tubes (post, courier service,...) | Storage of tubes upright (for approx. 1 hour) prior to further transport/dispatch |
| Recentrifugation of Sep tubes | Avoid absolutely! Could lead to a change in the analysis results (i.e. potassium) |