



User height
700–1700 mm



Max. load
77 kg



Environment
Indoor & outdoor

Rodeo

The Convaid Rodeo is a popular pediatric tilt-in-space wheelchair that is convenient, lightweight and folds compactly. The Rodeo provides up to 45° of adjustable tilt that offers pressure relief and upper body support. This tilt-in-space Rodeo is fully customizable to meet user's specific therapeutic needs and provides growth option through its adjustable seat depth feature.

Highlights about Rodeo



Adjustable tilt and back recline

Provides pressure relief and support.



Customisable growth capabilities

Adjustable seat depth extends chair life for a growing child.



Lightweight and durable

Lightweight frame combined with high Max. load - up to 77 kg.



Individual positioning

Accommodates functional or rest positions and increased sitting tolerance for maximum comfort.



Standard

5 to 45 degrees of adjustable tilt, recline, central mounted foot brake, adjustable one-piece push handle, closure strap, seat depth growth of 13 cm, frontwheels with directional stabilizer, 2-point positioning belt, elevating legrest /one-Piece, angle adjustable footplate.

Options

Size:	Brake	Tyres
- RD12	- Foot brake	- Solid
- RD14	- Attendant hand brakes	- Air
- RD16		- Solid w/transp. fitting
		- Air w/transp. fitting

Cushion colours



Blue



Red



Black



Pink



Purple



Green

Frame colour



Silver













Standard Rodeo RD14



Folded Rodeo

Measurements

		RD 12	RD 14	RD 16
	Effective seat width	305 mm	355 mm	405 mm
	Effective seat depth	230-360 mm	280-410 mm	330-460 mm
	Backrest height	560 mm	610 mm	635 mm
	Total mass	15 kg	16 kg	17 kg
	Max. load regular/transit	34/30 kg	46/46 kg	77/64 kg
	Overall width	580 mm	630 mm	680 mm
	Length	890 mm	890 mm	910 mm
	Overall height	1090 mm	1120 mm	1140 mm
	Folded length	860 mm	860 mm	1015 mm
	Folded width	585 mm	635 mm	685 mm
	Folded height	430 mm	460 mm	460 mm

Further measurements at our website

Nice to know

- Rodeo is suitable for walks in uneven terrain and nature.
- It is easy to adjust most of the settings and it is possible to use it for different children.



Accessories

1



Extended Headrest cover

With or without windows (Canopy)

12 C904745/ C904729
14 C904746/ C904730
16 C904747/ C904731

2



Rain cover and Mosquito net

Rain cover:

12/14 C904939 12/14 C904655
16 C904940 16 C904656

Mosquito net:

12/14 C904655
16 C904656

3



Headrest Extension

Required for transit model

12 C902863
14 C900484
16 C900138

4



Headwings

Padded:

12 C902673
14 C902653
16 C903147
Occi: C904779

5



H-harness with padded covers

With padded covers
C903988

6



3-pt positioning belt

With adjustable crotch strap
C904017

7



Full torso vest

12 C903261
14 C903260
16 C903259

8



Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5 cm smaller.

12 C904211
14 C904110

9



Swing away laterals - contoured

Small C904071 / C904466
Medium C904072 / C904467
Large C904073 / C904468

10



Cushions

Support, position and align cushion available.
12 C904350 / C904351 / C904352
14 C904353 / C904354 / C904355
16 C904356 / C904357 / C904358

11



Calf panel

C904361

12



Thigh support

Lateral: Medial:
12 C903668 12 C902668
14 C903648 14 C902647
16 C903141 16 C903140


13



Foot positioners

Small (18-25 cm) C903415
Medium (20-28 cm) C903416
Large (23-30 cm) C903417

14



Legrest extension

12 C904095
14 C904096
16 C904097


15



**Height adjustable
flip-up armrest**

C904346

16



Solid/Mesh side panels

12 C902672/ C902680
14 C902652/ C902681
16 C903145/ C903146

17



Rear anti-tip tubes (pair)

C904649

18



**Upper extremity support surface
(clear tray)**

Flip-up arm rest required
12 C903405
14 C903401
16 C903402


19



Upper extremity support tray

Foam tray
12 C903669
14/16 C903670

20



Anchor transit kit

C997124

21



Medical necessity storage

12 C904476
14/16 C904477

22



Under seat storage basket

12 C904529
14/16 C904530

23



Oxygen tank bag

C902870

24



Utility bag

C904531