

JOBST® FarrowWrap® 4000 and JOBST® FarrowHybrid®

Below knee A-D
Ready-to-wear order form




TO ORDER PLEASE EMAIL: info@hc21.ie

I consent to all personal data contained on this prescription form being shared with Healthcare 21 & Essity/ BSN Medical for the purpose of supply and delivery of this product. If applicable, I also consent to certain special categories of personal data, including but not limited to, my medical card number and relevant medical information being shared with HSE for funding purposes.

Signature: _____

Date: _____ Purchase Order No.: _____ Patient Name: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____



Legpiece - simple measuring points

1. Measure the calf and ankle circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing the oedema)
3. Measure the length of the leg posteriorly (B-D) and follow the contour of the leg. Select regular or tall

Left: Right:

Left: Right:

Left: Right:

†Measurement only required when ordering additional JOBST® FarrowHybrid®

When you order a JOBST® FarrowWrap® 4000 you will receive 1 x JOBST® FarrowHybrid® (20-30mmHg) wide

	X Small	Small	Medium	Large
JOBST FarrowWrap 4000 Legpiece				
cB Ankle	17-23 cm	19-28 cm	25-37 cm	29-43 cm
cC Calf	25-38 cm	32-48 cm	37-56 cm	43-67 cm
B-D Regular	30-33 cm	30-33 cm	30-33 cm	30-33 cm
B-D Tall	34-37 cm	34-37 cm	34-37 cm	34-37 cm

	X Small	Small	Medium	Large
JOBST® FarrowWrap® 4000				
JOBST® FarrowHybrid® size received with wrap	Medium-wide	Medium-wide	Large-wide	Large-wide*

*Do not apply JOBST® FarrowHybrid® if calf circumference exceeds 60 cm. When a large JOBST® FarrowWrap® 4000 is ordered, an additional non-compressive liner will be provided to accommodate calf circumference greater than 60 cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance.

When you order a JOBST® FarrowWrap® 4000 you will receive 1 x JOBST® FarrowHybrid® (20-30mmHg) wide

JOBST® FarrowWrap® 4000 Legpiece

Size: X Small Small Medium Large

Length: Regular* Tall**

Colour: Tan Black

Code: -

Quantity Required:

Repeat prescription required every months

*Regular: 30-33 cm **B-D** **Tall: 34-37 cm **B-D**

If an alternative JOBST® FarrowHybrid® size is required please order separately

JOBST® FarrowHybrid®

Shoe Size: Small (2-5) Medium (6-8) Large (9-11) X-Large (12-13)

Width: Regular* Wide**

Compression: 20-30mmHg - Taupe Black***
30-40mmHg - Sand

Code: -

Quantity Required:

Repeat prescription required every months

*Regular: 20-27 cm **A1** **Wide: 28-35 cm **A1** ***Available in wide variant only

JOBST® FarrowWrap® Order Codes

JOBST® FarrowWrap® 4000 Legpiece

Code	Size	Length	Colour
76661-00	Extra small	Regular	Beige
76661-08	Extra small	Regular	Black
76661-04	Extra small	Tall	Beige
76661-12	Extra small	Tall	Black
76661-01	Small	Regular	Beige
76661-09	Small	Regular	Black
76661-05	Small	Tall	Beige
76661-13	Small	Tall	Black
76661-02	Medium	Regular	Beige
76661-10	Medium	Regular	Black
76661-06	Medium	Tall	Beige
76661-14	Medium	Tall	Black
76661-03	Large	Regular	Beige
76661-11	Large	Regular	Black
76661-07	Large	Tall	Beige
76661-15	Large	Tall	Black



JOBST® FarrowHybrid®* - supplying compression to the foot and ankle

Code	Size	Colour
ADI (20-30mmHg)		
76664-00	Small (shoe size 2-5)	Beige
76664-01	Medium (shoe size 6-8)	Beige
76664-02	Large (shoe size 9-11)	Beige
76664-03	Extra large (shoe size 12-13)	Beige
76664-12	Small - wide (shoe size 2-5)	Beige
76664-13	Medium - wide (shoe size 6-8)	Beige
76664-14	Large - wide (shoe size 9-11)	Beige
76664-15	Extra large - wide (shoe size 12-13)	Beige
76664-16	Small - wide (shoe size 2-5)	Black
76664-17	Medium - wide (shoe size 6-8)	Black
76664-18	Large - wide (shoe size 9-11)	Black
76664-19	Extra large - wide (shoe size 12-13)	Black
ADII (30-40mmHg)		
76664-04	Small (shoe size 2-5)	Beige
76664-05	Medium (shoe size 6-8)	Beige
76664-06	Large (shoe size 9-11)	Beige
76664-07	Extra large (shoe size 12-13)	Beige
76664-08	Small - wide (shoe size 2-5)	Beige
76664-09	Medium - wide (shoe size 6-8)	Beige
76664-10	Large - wide (shoe size 9-11)	Beige
76664-11	Extra large - wide (shoe size 12-13)	Beige



*Do not apply JOBST® FarrowHybrid® if calf circumference exceeds 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance.