

JOBST® FarrowWrap®
A-D Below Knee

Ready-to-wear order form

TO ORDER PLEASE EMAIL: info@hc21.ie





I consent to all personal data contained on this prescription form being shared with Healthcare 21 & Essity/ BSN Medical for the purpose of supply and delivery of this product. If applicable, I also consent to certain special categories of personal data, including but not limited to, my medical card number and relevant medical information being shared with HSE for funding purposes.

Signature: _____

Date: _____ Purchase Order No.: _____ Patient Name: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____ Invoice Address: _____



Measuring instructions

Legpiece - 3 simple measuring points

1. Measure the widest calf **C** and narrowest ankle **B** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing the oedema)
3. Measure the length of the leg posteriorly from the floor to **D** measurement (2cm below the crease of the knee). Follow the contour of the leg. Use the sizing chart to select regular or tall to determine size

Footpiece - 2 simple measuring points

1. Measure **A¹** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column
3. Measure the **X** length of the foot. Use the sizing chart to select regular or long to determine size

	X Small	Small	Medium	Large	X Large
Legpiece					
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm
A-D Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm
NOTE: for A-D measure posteriorly and follow contours					
Footpiece					
A¹ Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
X Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm
X Tall	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm
NOTE: for X measure from the heel to the base of the little toe					

Note: 3cm allowance on either side of measurement.

☐ Legpiece

Size: ☐ X Small ☐ Small ☐ Medium ☐ Large ☐ X Large

Length A-D: ☐ Regular ☐ Tall

Fabric: ☐ Lite

☐ Classic

☐ Strong

BSN Code: -

Quantity Required:

Repeat prescription required every months

☐ Footpiece

Size: ☐ X Small ☐ Small ☐ Medium ☐ Large ☐ X Large

Length X: ☐ Regular ☐ Long

Fabric: ☐ Lite

☐ Classic

☐ Strong

BSN Code: -

Quantity Required:

Repeat prescription required every months

☐ JOBST FarrowHybrid

Size: ☐ Small ☐ Medium ☐ Large ☐ X Large

Width: ☐ Regular* ☐ Wide**

Compression:

20-30mmHg - ☐ Taupe

☐ Black***

BSN Code: -

Quantity Required:

Repeat prescription required every months

30-40mmHg - ☐ Sand

*Regular: 20-27 cm **A¹** **Wide: 28-35 cm **A¹** ***Available in wide variant only

A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap lower limb piece.

JOBST FarrowWrap Order Codes

JOBST FarrowWrap Lite (20-30mmHg)

Code	Size	Length	Colour
Legpiece			
76654-00	Extra Small	Regular	Beige
76654-05	Extra Small	Tall	Beige
76654-01	Small	Regular	Beige
76654-06	Small	Tall	Beige
76654-02	Medium	Regular	Beige
76654-07	Medium	Tall	Beige
76654-03	Large	Regular	Beige
76654-08	Large	Tall	Beige
76654-04	Extra Large	Regular	Beige
76654-09	Extra Large	Tall	Beige
Footpiece			
76652-05	Extra Small	Regular	Beige
76652-00	Extra Small	Long	Beige
76652-06	Small	Regular	Beige
76652-01	Small	Long	Beige
76652-02	Medium	Regular	Beige
76652-07	Medium	Long	Beige
76652-08	Large	Regular	Beige
76652-03	Large	Long	Beige
76652-09	Extra Large	Regular	Beige
76652-04	Extra Large	Long	Beige



JOBST FarrowWrap Classic (30-40mmHg)

Code	Size	Length	Colour
Legpiece			
76657-00	Extra Small	Regular	Beige
76657-05	Extra Small	Tall	Beige
76657-01	Small	Regular	Beige
76657-06	Small	Tall	Beige
76657-02	Medium	Regular	Beige
76657-07	Medium	Tall	Beige
76657-03	Large	Regular	Beige
76657-08	Large	Tall	Beige
76657-04	Extra Large	Regular	Beige
76657-09	Extra Large	Tall	Beige
Footpiece			
76656-05	Extra Small	Regular	Beige
76656-00	Extra Small	Long	Beige
76656-06	Small	Regular	Beige
76656-01	Small	Long	Beige
76656-07	Medium	Regular	Beige
76656-02	Medium	Long	Beige
76656-08	Large	Regular	Beige
76656-03	Large	Long	Beige
76656-09	Extra Large	Regular	Beige
76656-04	Extra Large	Long	Beige



JOBST FarrowWrap Strong (30-40mmHg)

Code	Size	Length	Colour
Legpiece			
76648-00	Extra Small	Regular	Beige
76648-05	Extra Small	Tall	Beige
76648-01	Small	Regular	Beige
76648-06	Small	Tall	Beige
76648-02	Medium	Regular	Beige
76648-07	Medium	Tall	Beige
76648-03	Large	Regular	Beige
76648-08	Large	Tall	Beige
76648-04	Extra Large	Regular	Beige
76648-09	Extra Large	Tall	Beige
Footpiece			
76646-04	Extra Small	Regular	Beige
76646-09	Extra Small	Long	Beige
76646-05	Small	Regular	Beige
76646-00	Small	Long	Beige
76646-06	Medium	Regular	Beige
76646-01	Medium	Long	Beige
76646-07	Large	Regular	Beige
76646-02	Large	Long	Beige
76646-08	Extra Large	Regular	Beige
76646-03	Extra Large	Long	Beige

JOBST FarrowHybrid* - supplying compression to the foot and ankle

Code	Size	Colour
ADI (20-30mmHg)		
76664-00	Small (shoe size 2-5)	Beige
76664-01	Medium (shoe size 6-8)	Beige
76664-02	Large (shoe size 9-11)	Beige
76664-03	Extra large (shoe size 12-13)	Beige
76664-12	Small - wide (shoe size 2-5)	Beige
76664-13	Medium - wide (shoe size 6-8)	Beige
76664-14	Large - wide (shoe size 9-11)	Beige
76664-15	Extra large - wide (shoe size 12-13)	Beige
76664-16	Small - wide (shoe size 2-5)	Black
76664-17	Medium - wide (shoe size 6-8)	Black
76664-18	Large - wide (shoe size 9-11)	Black
76664-19	Extra large - wide (shoe size 12-13)	Black

Code	Size	Colour
ADII (30-40mmHg)		
76664-04	Small (shoe size 2-5)	Beige
76664-05	Medium (shoe size 6-8)	Beige
76664-06	Large (shoe size 9-11)	Beige
76664-07	Extra large (shoe size 12-13)	Beige
76664-08	Small - wide (shoe size 2-5)	Beige
76664-09	Medium - wide (shoe size 6-8)	Beige
76664-10	Large - wide (shoe size 9-11)	Beige
76664-11	Extra large - wide (shoe size 12-13)	Beige
*Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm. Consider a JOBST FarrowWrap footpiece in this circumstance.		

