

Actimove[®]



Gauntlet

Wrist and Thumb Stabiliser


New universal size

Gauntlet

Wrist and Thumb Stabiliser

Secure MCP joint stabilisation, removable stay and universal fit for the left or right hand.

- Pain relief and healing are supported by stabilisation through the anatomically shaped stay
- MCP joint stabilisation through thumb stay while allowing full finger movement
- Great comfort from the soft **Latex-free** material and contoured strap ends
- Can be worn on the **left or right** hand
- Hand washable



Indications:

- Osteoarthritis of the thumb
- Irritation of the thumb
- Skier's thumb
- Sprains, strains
- Rheumatoid arthritis
- Post-operative care
- Post-traumatic care
- DeQuervain's tendosynovitis

Quick and easy to apply due to the slip-on sleeve and individually adjustable hook and loop fixings.



	Code	Size	Measurements	Colour
	72859-90	Small	13cm - 15cm	Black
	72859-91	Medium	15cm - 17.5cm	Black
	72859-92	Large	17.5cm - 20cm	Black
	72859-93	X Large	20cm - 23cm	Black