

Actimove®

Feel good. Move better.

NEW!

FOCUSED PROTECTION
AND STABILISATION
FOR THE LOWER LEG,
ANKLE AND FOOT

THE ACTIMOVE® WALKER RANGE



THERAPIES. HAND IN HAND.

BSN medical

HIGH QUALITY SUPPORT FOR EVERY STEP

ACTIMOVE® STANDARD WALKER, ACTIMOVE® STANDARD AIR WALKER AND ACTIMOVE® CLOSED SHELL AIR WALKER

The **Actimove Walker** range offers impressive features. The high quality, rigid outer skeleton ensures high protection of the foot, ankle and lower leg. Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient. Walking with **Actimove Walker** products feels natural through the rocker bottom which simulates a physiological gait.

••• Actimove Standard Walker

Secure and lightweight stabilisation for the lower leg, ankle and foot



••• Actimove Standard Air Walker

Adjustable stabilisation with integrated pump and valve system to easily alter desired compression level



••• Actimove Closed Shell Air Walker

Designed to offer the highest protection for the lower leg, ankle and foot

Indications

- Fractures of the foot and ankle
- Lower extremity soft-tissue injuries
- Severe ankle and foot sprains
- Post-operative and rehabilitation use

For Closed Shell Air Walker only:

- Achilles tendon ruptures (only in combination with separately available Actimove Heel Wedge for Walkers)

ADVANCED FEATURES FOR SAFETY AND PATIENT CARE



..... Circumferential strapping
For critical stability



..... Breathable liner
For high wearing comfort



..... Shock absorbing sole and
rocker bottom
Minimises impact and
simulates a natural gait

INTEGRATED PUMP SYSTEM OF ACTIMOVE STANDARD AIR WALKER AND ACTIMOVE CLOSED SHELL AIR WALKER



..... Integrated pump and
valve system
Designed to easily adjust
desired compression level



ADDITIONAL ADJUSTABLE FEATURES OF THE ACTIMOVE CLOSED SHELL AIR WALKER



..... Heel Wedge
Adjustable heel wedge for
Achilles tendon injuries



..... Adjustable toe guard
Provides protection in the
toe area

Ordering information

Walker		Size	BSN Code - Low	BSN Code - High
Actimove Standard Walker Low / High Lightweight walker		XS	76270-30	76270-35
		S	76270-31	76270-36
		M	76270-32	76270-37
		L	76270-33	76270-38
		XL	76270-34	76270-39
Actimove Standard Air Walker Low / High Lightweight walker with integrated pump and valve system		XS	76272-30	76272-35
		S	76272-31	76272-36
		M	76272-32	76272-37
		L	76272-33	76272-38
		XL	76272-34	76272-39
Actimove Closed Shell Air Walker Low / High Lightweight closed shell walker with integrated pump and valve system		XS	76276-10	76276-15
		S	76276-11	76276-16
		M	76276-12	76276-17
		L	76276-13	76276-18
		XL	76276-14	76276-19
Walker Accessories		Size	BSN Code	
Actimove Walker Heel Wedge Heel wedge for Achilles tendon injuries for use with the Closed Shell Air Walker		Universal	76277-94	



Size Guide

Size	
XS	< 3
S	3.5 - 5.5
M	6 - 8.5
L	9 - 11.5
XL	> 12