

The Affix drop foot AFO

Affix Ankle Brace

The Affix is a lightweight ankle orthosis, designed with the active support for drop foot and similar conditions in mind.

Through the provision of immediate support as soon as the foot is lifted this simple brace provides a discernible improvement to the patient's gait. The uncomplicated nature of the design is also replicated in the Affix's simplicity to fit. Made from breathable fabric, the ankle wrap provides increased comfort when wearing for long periods of time, whilst also preventing irritation or hampering freedom of movement.



Features:

- Effective drop foot splint device in dual parts, supporting above the ankle and with a non-slip section fitting between the tongue and shoelaces of footwear.
- Fully adjustable elastic strap provides tailored support between dual sections.
- Targeted, localised gel padding ensures a discreet and comfortable fit.
- Lightweight and intuitive to fit.
- Shoeless wrap available for use without footwear (separate Affix Ankle Brace sti 11 required)
- Replacement traction elastic, S hooks and shoe inserts available separately.

Indications:

- Drop Foot
- Peroneal Neuropathies
- Dysfunctional Gait

Contraindications:

- Spastic Paralysis
- Circulatory/sensory disorders



NOTE:
Measure ankle circumference just above the ankle bone.

Affix Ankle Brace

SIZE	ANKLE CIRC.		AFFIX ANKLE BRACE		AFFIX CALF SUPPORT	
	INCHES	CM	BLACK	BEIGE	BLACK	BEIGE
X-Small	5¾ - 6¾	13 - 17	ABB-XS	ABBE-XS	ABB-CS-XS	ABBE-CS-XS
Small	6¾ - 8¼	17 - 21	ABB-S	ABBE-S	ABB-CS-S/M	ABBE-CS-S/M
Medium	8¼ - 10	21 - 25	ABB-M	ABBE-M	ABB-CS-S/M	ABBE-CS-S/M
Large	10 - 11½	25 - 29	ABB-L	ABBE-L	ABB-CS-L	ABBE-CS-L

Shoeless Wrap

SIZE	FOOT CIRC.		SHOELESS WRAP - LEFT		SHOELESS WRAP - RIGHT	
	INCHES	CM	BLACK	BEIGE	BLACK	BEIGE
Small	6¾ - 8¼	17 - 21	ABB-SW-S-L	ABBE-SW-S-L	ABB-SW-S-R	ABBE-SW-S-R
Medium	8¼ - 10	21 - 25	ABB-SW-M-L	ABBE-SW-M-L	ABB-SW-M-R	ABBE-SW-M-R
Large	10 - 11½	25 - 29	ABB-SW-L-L	ABBE-SW-L-L	ABB-SW-L-R	ABBE-SW-L-R

Replacement Accessories

SIZE	TRACTION ELASTIC		SHOE INSERTS (2 PK)		S HOOKS (2 PK)
	BLACK	BEIGE	BLACK	BEIGE	
Small	ABB-E-S	ABBE-E-S	ABB-SI	ABBE-SI	ABB-SH
Medium	ABB-E-M	ABBE-E-M			
Large	ABB-E-L	ABBE-E-L			

