



SCISSORS

CAST ROOM SCISSORS



MP1994 Supercut Cast Scissor

Lower blade with fine serrations and sharp- angled upper blade



MP1991C 19cm Cast Cutter - Long Blade Guarded End
MP1993C 23cm Cast Cutter - Long Blade Guarded End

Serrated bottom blade designed for specifically cutting synthetic casts.



MP2011 Utility scissors - standard

Plastic handles and plastic tips



MP1979 14cm Lister bandage scissor

MP1980 18cm Lister bandage scissor

MP1981 20cm Lister bandage scissor

High quality stainless steel



MP1992C 18cm Cast Cutter - Short Blade

Serrated bottom blade designed for specifically cutting synthetic casts.



MP1989-19C 19cm Eccentric Cast Cutter

MP1989-16C 16cm Eccentric Cast Cutter

Very strong cutting action and high quality stainless steel



MP1980TC 18cm Lister bandage scissor

Tungsten insert blades - High quality S/Steel



MP1984 24cm Bruns Bandage Scissors

Plain bottom blade

MP1985 24cm Bruns Bandage Scissors

Bottom blade serrated, very strong heavy duty cast cutter



MP1811TC 12cm Curved Nibbler with Tungsten insert blade

Tungsten insert blade bottom blade serrated



MP1811 12cm Curved Nibbler Scissor

Popular design with serrated bottom jaw



MP2022 - 18cm Lister Bandage scissor - Left hand

High quality stainless steel



MP2020 pack of 6

MP2020S Individual

Economy Plaster Scissor 18cm - quality scissor with serrated bottom blade - packs of 6

INSTRUMENTS



MP2010 Cast Spreader

Tri Blade cast spreader 23cm – very easy to use with fine tips to jaws.



MP2021 18cm Esmarch Plaster Knife

Stainless steel - one edge sharp



MP2034 28cm Henning Cast Spreader

Traditional pattern cast spreader



MP1996 23cm Stille Cast Shears

MP1997 26cm Stille Cast Shears

MP1998 37cm Stille Cast Shears

High quality precision cut



MP2032 18cm Wolff Cast Spreader

MP2033 24cm Wolff Cast Spreader

Deep cut serrations for a better grip



MP2016 - Ringcutter

Strong handles with hardened cutting blade



MP2035 18cm Charnley Wire Holders

Serrated jaws for removing wires

MP2036 18cm Charnley Wire Holders

Tungsten serrated jaws for removing wires

ABOUT MEDEZINE

For over 20 years Medezine have designed and manufactured high quality medical equipment, tailored to the clinical demands of the modern N.H.S.

Medezine successfully exports worldwide. Based in Sheffield, England, Medezine achieves and maintains high standards, having BS EN ISO 9001 accreditation. All products are CE marked in accordance with European directive 93/42/EEC.

Medezine offer a comprehensive range of casting equipment, which includes, high quality saws, scissors/instruments, and specialist trolleys.

Customer care is paramount throughout the whole organisation.

Medezine Limited,
Churchill Way,
Chapelton, Sheffield,
England, S35 2PY

Tel: +44 (0) 114 257 6002
Fax: +44 (0) 114 246 6607
Email: sales@medezine.co.uk
Web: www.medezine.co.uk