

Humeral Fracture Brace

Features:

- Lightweight design, constructed from trimmable polyethylene
- Allows full ROM at shoulder and elbow for early physiotherapy
- Lined with closed-cell foam for increased comfort
- Lining does not absorb moisture
- Adjustable hook and loop straps to maintain good compression
- Universal - fits right or left arm

Indications:

- Non-operative closed humeral fractures
- Humeral diaphyseal fractures
- Post external/internal fixation
- Deep humeral bone bruises/post biopsy

Size	Mid Humeral Circumference		Code
	Inches	cm	
X-Small	5 - 7½	12.7 - 19	HFB-XS
Small	7½ - 9½	19 - 24	HFB-S
Medium	9½ - 11½	24 - 29	HFB-M
Large	11½ - 13½	29 - 34	HFB-L
X-Large	13½ - 15½	34 - 39	HFB-XL



Over the Shoulder Humeral Fracture Brace

Features:

- Lightweight design, constructed from trimmable polyethylene
- Allows full ROM at elbow for early physiotherapy
- Proximal shoulder cap for increased stability
- Lined with closed-cell foam for increased comfort
- Lining does not absorb moisture
- Adjustable hook and loop straps to maintain good compression
- Universal - fits right or left arm

Indications:

- Higher fractures where shoulder ROM needs to be restricted
- Non-operative closed humeral fractures
- Humeral diaphyseal fractures
- Deep humeral bone bruises/post biopsy
- Post external/internal fixation



Size	Mid Humeral Circumference		Code
	Inches	cm	
X-Small	5 - 7½	12.7 - 19	HFB-OS-XS
Small	7½ - 9½	19 - 24	HFB-OS-S
Medium	9½ - 11½	24 - 29	HFB-OS-M
Large	11½ - 13½	29 - 34	HFB-OS-L
X-Large	13½ - 15½	34 - 39	HFB-OS-XL