

Shell Walker/ Shell Air Walker

Features:

- Substitutes ankle motion and aids propulsion whilst decreasing heel-strike force on calcaneus
- Includes inflatable aircells for further support while reducing swelling and evenly distributing pressure (**Shell Air Walker only**)
- Contoured strut design provides better fit and support
- Rocker bottom assists gait
- Front and back panels for superior protection
- Available in standard (17") or short (11") height
- Replacement liners available separately

Indications:

- Grade II/III ankle sprains
- Stable fractures of the fibia, talus, calcaneus and medial malleolus
- Post ligament, soft tissue and tendon surgery
- Forefoot or midfoot injury



Shell Air Walker Pump



Short Version



Size	Shell Walker						Shell Air Walker			
	Shoe Size		Height		Replacement Liners		Height		Replacement Liners	
	UK	Euro	Standard	Short	Standard	Short	Standard	Short	Standard	Short
Small	< 5	< 38	SW-S	SWS-S	SW-LIN-S	SWS-LIN-S	SW-A-S	SWS-A-S	SW-A-LIN-S	SWS-A-LIN-S
Medium	6 - 8	38 - 42	SW-M	SWS-M	SW-LIN-M	SWS-LIN-M	SW-A-M	SWS-A-M	SW-A-LIN-M	SWS-A-LIN-M
Large	9 - 11	43 - 46	SW-L	SWS-L	SW-LIN-L	SWS-LIN-L	SW-A-L	SWS-A-L	SW-A-LIN-L	SWS-A-LIN-L
X-Large	12 - 14	47 - 50	SW-XL	SWS-XL	SW-LIN-XL	SWS-LIN-XL	SW-A-XL	SWS-A-XL	SW-A-LIN-XL	SWS-A-LIN-XL