




COMVI™

Enteral stent (Flare)



For malignant colorectal / gastric outlet obstruction


Taewoong Niti-S™ Comvi™ Enteral Stent (Flare)

for malignant colorectal / gastric outlet obstruction



Features

- **Triple layered construction** : Biocompatible PTFE membrane tube is held between inner and outer mesh
- Uncovered proximal flare end to reduce migration
- Large cell structure of body with thick wire to reduce the risk of fracture and enhance radial force
- PTFE membrane prevents the risk of the tissue invasion
- **Large diameter (up to 26mm)** is loaded into 10.5Fr delivery system
- **Minimum foreshortening** for accurate stent placement



- **Radiopaque marker** : 3(three) at Flared end and 2(two) on both ends of body

Ordering Information

TTS (Through The Scope)

► Enteral Colonic stent

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CCT1806P	18	6	10.5	220
CCT1808P		8		
CCT1810P		10		
CCT1812P		12		
CCT2006P	20	6		
CCT2008P		8		
CCT2010P		10		
CCT2012P		12		
CCT2206P	22	6		
CCT2208P		8		
CCT2210P		10		
CCT2212P		12		
CCT2406P	24	6		
CCT2408P		8		
CCT2410P		10		
CCT2412P		12		
CCT2606P	26	6		
CCT2608P		8		
CCT2610P		10		
CCT2612P		12		

► Pyloric / duodenal stent

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806P	18	6	10.5	180
PCT1808P		8		
PCT1810P		10		
PCT1812P		12		
PCT2006P	20	6		
PCT2008P		8		
PCT2010P		10		
PCT2012P		12		
PCT2206P	22	6		
PCT2208P		8		
PCT2210P		10		
PCT2212P		12		



Aquilant Interventional

UK Distributors

Aquilant Interventional
Aquilant House, Units B1-B2
Bond Close, Kingsland Business Park
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007
F: +44 (0) 1256 365486
E: contactus@aquilantinterventional.com
W: www.aquilantinterventional.com

