





COMVI™ Enteral colonic stent

for malignant colorectal obstruction



Feature

- **Triple layered construction** : Biocompatible PTFE membrane tube is held between inner and outer mesh
 - Unfixed cell structure enables stent to conform to tortuous anatomy
 - PTFE membrane prevents the risk of the tissue invasion
 - Outer wire mesh prevents the risk of migration
- **Minimum foreshortening** for accurate stent placement
- **Radiopaque marker** : 4(four) at both covered part ends



Ordering Information

TTS (Through The Scope)								
Code	Stent		Delivery					
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)				
CCT1806BA	18	6	10.5	220				
CCT1808BA		8						
CCT1810BA		10						
CCT1812BA		12						
CCT2006BA	20	6			10.5	220		
CCT2008BA		8						
CCT2010BA		10						
CCT2012BA		12						
CCT2206BA	22	6					10.5	220
CCT2208BA		8						
CCT2210BA		10						
CCT2212BA		12						

OTW (Over The Wire)								
Code	Stent		Delivery					
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)				
CC1806BA	18	6	14	70				
CC1808BA		8						
CC1810BA		10						
CC1812BA		12						
CC2006BA	20	6			14	70		
CC2008BA		8						
CC2010BA		10						
CC2012BA		12						
CC2206BA	22	6					16	70
CC2208BA		8						
CC2210BA		10						
CC2212BA		12						
CC2406BA	24	6	16	70				
CC2408BA		8						
CC2410BA		10						
CC2412BA		12						
CC2606BA	26	6			18	70		
CC2608BA		8						
CC2610BA		10						
CC2612BA		12						
CC2806BA	28	6					18	70
CC2808BA		8						
CC2810BA		10						
CC2812BA		12						

Released Articles

* **Comparison of efficacies between stents for malignant colorectal obstruction: a randomized, prospective study**

GASTROINTESTINAL ENDOSCOPY, Vol. 72, No. 2: 2010 Semi Park

* **Comparison of a Newly Designed Double-Layered Combination Covered Stent and D-Weave Uncovered Stent for Decompression of Obstructive Colorectal Cancer: A Prospective Multicenter Study** The ASCRS Vol. 53, No. 8: 2010 Chang Mo Moon, MD, et al.



Aquilant Interventional

UK Distributors

Aquilant Interventional
Aquilant House, Units B1-B2
Bond Close, Kingsland Business Park
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007
F: +44 (0) 1256 365486
E: contactus@aquilantinterventional.com
W: www.aquilantinterventional.com