






COMVI™ Pyloric / duodenal stent

for malignant gastric outlet obstruction



Feature

- **Triple layered construction** : Biocompatible PTFE membrane tube is held between inner and outer mesh
 - Unfixed cell structure enables stent to conform to tortuous anatomy
 - PTFE membrane prevents the risk of tissue invasion
 - Outer wire mesh prevents the risk of migration
- **Minimum foreshortening** for accurate stent placement
- **Radiopaque marker** : 4(four) at both covered part ends



◀ After the endoscope was withdrawn, the Comvi™ stent showed good conformability along the duodenum

[GASTROINTESTINAL ENDOSCOPY Vol. 75, Issue. 4, p757-763]

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806BA	18	6	10.5	180
PCT1808BA		8		
PCT1810BA		10		
PCT1812BA		12		
PCT2006BA	20	6		
PCT2008BA		8		
PCT2010BA		10		
PCT2012BA		12		
PCT2206BA	22	6		
PCT2208BA		8		
PCT2210BA		10		
PCT2212BA		12		

* Available 220cm delivery system

Coding : PCT**__BA-22

OTW (Over The Wire)						
Code	Stent		Delivery			
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		
PC1806BA	18	6	14	135		
PC1808BA		8				
PC1810BA		10				
PC1812BA		12				
PC2006BA	20	6				
PC2008BA		8				
PC2010BA		10				
PC2012BA		12				
PC2206BA	22	6			16	135
PC2208BA		8				
PC2210BA		10				
PC2212BA		12				
PC2406BA	24	6				
PC2408BA		8				
PC2410BA		10				
PC2412BA		12				
PC2606BA	26	6	18	135		
PC2608BA		8				
PC2610BA		10				
PC2612BA		12				
PC2806BA	28	6				
PC2808BA		8				
PC2810BA		10				
PC2812BA		12				

Released Articles

* **Management of malignant gastric outlet obstruction with a modified triple-layer covered metal stent**

GASTROINTESTINAL ENDOSCOPY 2012 Hiroyuki Isayama, MD, PhD, et al.

* **Self-expanding metallic stents for gastric outlet obstruction resulting from stomach cancer: a preliminary study with a newly designed double-layered pyloric stent** GASTROINTESTINAL ENDOSCOPY Vol. 66, No.6: 2007 Sum Mi Lee, MD, PhD, et al.



Aquilant Interventional

UK Distributors

Aquilant Interventional
Aquilant House, Units B1-B2
Bond Close, Kingsland Business Park
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007
F: +44 (0) 1256 365486
E: contactus@aquilantinterventional.com
W: www.aquilantinterventional.com