




# **CONIO™** Esophageal stent



# Taewoong Niti-S™ Conio™ Esophageal stent for hypopharyngeal strictures

## Features

- **Small diameter with proximal head design** specially designed for refractory hypopharyngeal strictures
- **Silicone covering** : Reduce the risk of tumor ingrowth
- **Visible green suture for easy removal**
- **Radiopaque marker** : 4(four) at both ends & 2(two) in the middle
- ✳ **Dr. Massimo Conio invented** this stent and has treated patients with refractory hypopharyngeal strictures after combined therapy for laryngeal cancer

## Ordering Information

Distal Release					Proximal Release					
Code	Stent		Delivery		Code	Stent		Delivery		
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)	
ES1006FN	10	6	12	70	ESP1006FN	10	6	14	70	
ES1008FN		8			ESP1008FN		8			
ES1010FN		10			ESP1010FN		10			
ES1012FN		12			ESP1012FN		12			
ES1014FN		14			ESP1014FN		14			
ES1015FN		15			ESP1015FN		15			
ES1206FN	12	6	14		ESP1206FN	12	6			14
ES1208FN		8			ESP1208FN		8			
ES1210FN		10			ESP1210FN		10			
ES1212FN		12			ESP1212FN		12			
ES1214FN		14			ESP1214FN		14			
ES1215FN		15			ESP1215FN		15			
ES1406FN	14	6	14		ESP1406FN	14	6	16		
ES1408FN		8			ESP1408FN		8			
ES1410FN		10			ESP1410FN		10			
ES1412FN		12			ESP1412FN		12			
ES1414FN		14			ESP1414FN		14			
ES1415FN		15			ESP1415FN		15			
ES1606FN	16	6	20		ESP1606FN	16	6			16
ES1608FN		8			ESP1608FN		8			
ES1610FN		10			ESP1610FN		10			
ES1612FN		12			ESP1612FN		12			
ES1614FN		14			ESP1614FN		14			
ES1615FN		15			ESP1615FN		15			

## Released Articles

- \* **Temporary placement of a fully covered self expanding metal stent to allow therapeutic ERCP**  
Endoscopy 2014; 46: E419 Mangiavillano Benedetto et al.
- \* **Self-Expanding Stents in Benign Esophageal Strictures**  
Gastrointestinal Endoscopy doi:10.1016/j.tgie/2008.07.001 Massimo Conio, MD, et al.
- \* **A modified self-expanding Niti-S stent for the management of benign hypopharyngeal strictures**  
GASTROINTESTINAL ENDOSCOPY Vol. 65, No. 4: 2007 Massimo Conio



Aquilant Interventional

### UK Distributors

Aquilant Interventional  
Aquilant House, Units B1-B2  
Bond Close, Kingsland Business Park  
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007  
F: +44 (0) 1256 365486  
E: [contactus@aquilantinterventional.com](mailto:contactus@aquilantinterventional.com)  
W: [www.aquilantinterventional.com](http://www.aquilantinterventional.com)

