



The **Dale Abdominal Binder**, specifically designed for post-operative use, is made entirely of Dale's exclusive stretch material, which permits hook and loop closure at any point for evenly distributed compression.

Dale®

Always Better

pat. 8,066,657
CE
ABDOMINAL BINDER
Dale
418
60 - 75 IN / 152 - 191 CM
1-508-695-9316 • MADE IN U.S.A.
Not made with natural rubber latex
REV 1



Abdominal Binder with EasyGrip™ Strip

THE STRETCH COMPRESSION ALLOWS FOR SUDDEN LUNG EXPANSION, WHICH IS NECESSARY FOR COUGHING AND DEEP BREATHING THERAPY.

Studies¹ show that a post-operative regimen that includes the Dale Abdominal Binder with EasyGrip Strip encourages better pulmonary function, allows earlier ambulation, and alleviates incisional pain.

The Dale Abdominal Binder with EasyGrip Strip splints the wound and encourages patients to breathe deeply, turn in bed, cough to clear their lungs, and ambulate immediately following surgery. As a result, there is less risk of sedentary complications such as atelectasis, pneumonia, and phlebitis.



- EasyGrip Strip holds drainage bulbs in place.
- Can be cut to allow for tubing and will not fray or unravel.
- Made of exclusive stretch material with hook that can be engaged at any point along the binder.
- Innovative paneled construction helps prevent riding, roping, and folding over
- Evenly distributes compression.
- Facilitates deep breathing.
- Not made with natural rubber latex.



APPLICATIONS

The Dale Abdominal Binder with EasyGrip Strip provides support and comfort following procedures such as:

- Cesarean Section
- Gastroplasty
- Hysterectomy
- Abdominoplasty
- Liposuction
- Childbirth

ORDERING INFORMATION		
Product Number	Size	Quantity
525	6" Wide, Stretches to fit 72cm-127cm (28"-50")	1 per box
410	9" Wide, Stretches to fit 76cm-114cm (30"-45")	1 per box
411	9" Wide, Stretches to fit 117cm-157cm (46"-62")	1 per box
418	9" Wide, Stretches to fit 152cm-191cm (60"-75")	1 per box
810	12" Wide, Stretches to fit 76cm-114cm (30"-45")	1 per box
811	12" Wide, Stretches to fit 117cm-157cm (46"-62")	1 per box
818	12" Wide, Stretches to fit 152cm-191cm (60"-75")	1 per box
820	12" Wide, Stretches to fit 183cm-239cm (72"-94")	1 per box
920	15" Wide, Stretches to fit 183cm-239cm (72"-94")	1 per box

HCPCS Code A4649

Dale
www.dalemed.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
EasyGrip™ is a trademark of Dale Medical Products, Inc. in the USA.
©2017 Dale Medical Products, Inc.

¹Kuan TA, Serrett C, and Ali J; The Effect of Abdominal Binder of Postoperative Pulmonary Function. Infections in Surgery 1983;2(11):875-883