

PATIENT COMFORT AND  
SATISFACTION SHOULD  
NEVER BE COMPROMISED

# think

of the ramifications

## DALE® POST-SURGICAL BRA

### CLINICAL ADVANTAGES

Features	Benefits	Clinical Advantages
<b>Seamless, One-Piece Design</b>	<ul style="list-style-type: none"><li>• Provides maximum comfort throughout recovery</li><li>• Soft, stretchy material conforms to breasts</li><li>• Natural looking shape under clothing</li><li>• Easy to apply to patient or self</li><li>• Comfortable straps prevent shoulder irritation or "digging"</li></ul>	<ul style="list-style-type: none"><li>• Secures dressing without the need for tape or "ace-type" bandages</li><li>• Smooth, seamless design prevents additional trauma to surgical site</li><li>• Provides needed compression and support.</li></ul>
<b>Exclusive Stretch Material</b>	<ul style="list-style-type: none"><li>• Secure, soft comfort for duration of recovery</li><li>• Lightweight material</li><li>• Machine washable</li></ul>	<ul style="list-style-type: none"><li>• Stretch knit fabric permits surgical site to breathe and won't cause skin irritation</li><li>• Material expands during recovery period to accommodate edema</li></ul>
<b>Quick-Fastening Hook &amp; Loop Front Closure</b>	<ul style="list-style-type: none"><li>• Easy application</li><li>• No cumbersome hook and eye closure or clips.</li></ul>	<ul style="list-style-type: none"><li>• Facilitates dressing changes</li><li>• Easy examination of chest</li><li>• Easy access to dressings</li><li>• Closing the bra in the front rather than in the back may reduce arm discomfort.</li></ul>
<b>Convenient Sizing</b>	<ul style="list-style-type: none"><li>• Reduces inventory and cost with four sizes available: M, L, XL, XXL</li><li>• Only two bras required per patient; one can be washed while wearing the other</li></ul>	<ul style="list-style-type: none"><li>• Proper sizing ensures full coverage and gentle support</li></ul>
<b>Improves Patient Comfort</b>	<ul style="list-style-type: none"><li>• Front closure</li><li>• Soft, stretch fabric</li><li>• Can be worn as a leisure or sleep bra following recovery.</li></ul>	<ul style="list-style-type: none"><li>• Enhances recovery by increasing patient confidence and sense of security</li><li>• Encourages normal daily activities with increased support and comfort</li><li>• Transitional bra between surgery and regular bra</li></ul>
<b>Not Made With Natural Rubber Latex</b>	<ul style="list-style-type: none"><li>• Eliminates patient and care-giver latex reactions</li></ul>	<ul style="list-style-type: none"><li>• Prevents latex hypersensitivity</li></ul>



800-343-3980  
[www.dalemed.com](http://www.dalemed.com)

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.  
©2014 Dale Medical Products, Inc.  
CA-BRA\_Rev\_C