



**Dale Tracheostomy Tube Holders** are designed to provide secure positioning and to minimize movement of the tracheostomy tube for all your patients.



## Tracheostomy Tube Holders

WHETHER YOU CHOOSE THE WELL-KNOWN DALE 240 BLUE® FOR ADULTS, PEDIPRINTS® SERIES FOR KIDS, OR THE BARIATRIC 244 BLUE® YOU CAN BE SURE IT'S THE RIGHT CHOICE.

Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheoesophageal fistula, tracheal stenosis, or airway granuloma. That's why it is important to know that no other securement method stabilizes tracheostomy tubes better than Dale. Now you can eliminate the frustrations associated with twill ties and other holders while minimizing secondary complications at the same time.

- Secure fastener tabs limit movement of the tracheostomy tube reducing accidental decannulation.
- Exclusive moisture wicking neckband lining reduces the risk of skin breakdown by minimizing the amount of moisture absorbed.
- Adjustable two piece neckband.
- PediPrints® series offers kids fun, comforting prints.†
- Fastener tabs quickly and easily attach to any size tracheostomy tube.
- Stretch material‡ allows for cough reflex, accommodates edema, and ensures a snug fit without compromising security.
- Not made with natural rubber latex.

† On the 241 and 242 tube holder  
‡ On the 240 Blue® tube holder

ORDERING INFORMATION		
Product Number	Size	Quantity
240	240 Blue® 2.54cm(1") wide, fits up to 49.5cm(19.5") neck	10 per box
241	241 PediPrints® 2.54cm(1") wide, fits up to 45.5cm(18") neck	10 per box
242	242 PediPrints® 1.91cm(0.75") wide, fits up to 23cm(9") neck	10 per box
244	244 Blue® 3.81cm(1.5") wide, fits up to 71cm(28") neck	10 per box
200	Dale One Piece Holder 2.54cm(1") wide neckband, fits 24.1cm - 45.5cm(9.5" - 18") neck	30 per box

HCPCS Code A7526

**Dale**  
www.dalemed.com

Dale®, 240 Blue®, 244 Blue®, & PediPrints® are registered trademarks of Dale Medical Products, Inc. in the USA.  
©2017 Dale Medical Products, Inc.