

Elbur Residential Care Bed


- ▶ Modern and homely look with all the key functions of an acute profiling bed
- ▶ Split side rails and central braking systems available
- ▶ Choice of woods available
- ▶ Low height of 26cm standard
- ▶ Single and extra wide widths available
- ▶ Designed and manufactured in Europe
- ▶ 5 year warranty



General Information:

Save work load:	175 kg
Max. User weight:	135 kg
Total weight:	124,1 kg
Heaviest single part:	30,2 kg
Bed Dimensions:	Single Length: 207cm, Width: 90cm, Extra width: 105cm
Back part length:	80cm
Heading adjustability:	70°
Foot adjustability:	20°
Height under mattress support:	13cm
Caster size:	100 mm always with brake.

Colours:



R5111
Beech

R3114
Nabucco

H3704
Walnut

R4272
Oak

R4991
Cherry

R5052
Acacia

R5311
Light Beech

R5703
Maple

For further information please contact

Greg Ward - Nursing Care | Business Unit Manager on **086 815 9372**

Stephen Dennehy - Key Account Manager on **087 678 1287**

or info@hc21.ie