


Eleview™

LIFT THAT LASTS

>45

**MINUTES
LIFTING**



**LONG LASTING
PERFORMANCE
EXCELLENT
VISIBILITY**

Eleview is a trademark of Cosmo Technologies Ltd.
Distributed by **Aquilant Ltd**

SUBMUCOSAL INJECTABLE COMPOSITION – DESIGNED FOR EASY AND SAFE RESECTION PROCEDURES

INTENDED USE

Elevview™ is a new submucosal injectable composition intended for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early stage cancers, or other gastrointestinal mucosal lesions, prior to excision with a snare or endoscopic device.



5 ampoules / pouch (10 ml / ampoule)

PERFORMANCE



- Long-lasting cushion with > 45 min lifting time
- Appropriate for challenging polyps regardless of size, location, or type
- Improved margin visualisation for less potential risk of perforation

EFFICIENCY



- Requires less volume to create cushions
- Fewer reinjections and piecemeal excisions compared to saline that can increase your time saving

VISIBILITY



- Dyed with Methylene Blue for improved visualisation of target lesion margins

SAFE AND EASY



- Pre-mixed, sterile and ready to use
- Low-viscosity emulsion fits through a 23 G needle and is easy to inject

KEY FEATURES



Aquilant

Customer Support Enquiries Telephone: +44 (0)1256 365 456
contactus@aquilantendoscopy.com
www.aquilantendoscopy.com

Elevview is a trademark of Cosmo Technologies Ltd,
Riverside II – Sir John Rogerson's Quay, Dublin 2, IE
info@cosmotechnologies.com
Manufacturer: Cosmo Technologies Ltd., IE

Distributed by **Aquilant Ltd**