



Product Catalog

October 2021

www.reachsurgical.com



CONTENTS

Endoscopic Stapling 4-7

Open Stapling 8-9

Energy 10-11

EXPAND

THE SCOPE OF SURGERY

Reach Surgical, a Genesis Medtech company, is committed to serving high-quality products to clinicians and patients worldwide.

We deliver advanced technology and cost-effective minimally invasive surgical products to medical institutions and provide patients with reliable and economical treatment solutions.

ENDO III Stapler



Series	Product Code	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO III 45mm Articulating Stapler	REC45A S	282	12	3
	REC45A	332		
	REC45A L	432		
ENDO III 60mm Articulating Stapler	REC60A S	298		
	REC60A	348		
	REC60A L	448		


iReach Magnum Powered Stapler

Product Code	Staple Line Length (mm)	Shaft Length (mm)	Instrument Length (mm)	Trocar Compatibility	Max number of reloads	QTY/BOX
IM60AS	60	298	546	12 mm	12	3
IM45AS	45	282	530			
IM60AM	60	348	596			
IM45AM	45	332	580			
IM60AL	60	448	696			
IM45AL	45	432	680			

ENDO III Cartridge



Product Code	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Number of staple rows	QTY/BOX
REC 45 GRA	Gray	2.0	0.75	Extra Thin	6	12
REC 60 GRA						
REC 45 WHT	White	2.5	1.0	Vascular		
REC 60 WHT						
REC 45 BLU	Blue	3.5	1.5	Regular		
REC 60 BLU						
REC 45 GLD	Gold	3.8	1.75	Regular/Thick		
REC 60 GLD						
REC 45 GRN	Green	4.1	2.0	Thick		
REC 60 GRN						
REC 45 BLK	Black	4.4	2.2	Extra Thick		
REC 60 BLK						



ENDO SRC Universal Stapler



Product Code	Shaft Length (mm)	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO SRC S	65	696	25	3
ENDO SRC	155	680		
ENDO SRC L	255	596		

ENDO SRC Single Use Loading Unit



Product Code *	Staple Line Length (mm)	Cartridge Color	Convex Points	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Trocar Compatibility	Number of staple rows	QTY/BOX
ENDO SRC4525 L ENDO SRC4525 R / ENDO SRC4525 BR	45	White	N/A	2.5	1.0	Thin	12 mm	6	6
ENDO SRC6025 L ENDO SRC6025 R / ENDO SRC6025 BR	60								
ENDO SRC4535 L ENDO SRC4535 R / ENDO SRC4535 BR	45	Blue	Yes	3.5	1.5	Regular			
ENDO SRC6035 L ENDO SRC6035 R / ENDO SRC6035 BR	60								
ENDO SRC6040 L ENDO SRC6040 R / ENDO SRC6040 BR	60	Gold	Yes	4.0	1.75	Regular/ Thick			
ENDO SRC4548 L ENDO SRC4548 R / ENDO SRC4548 BR	45	Green	Yes	4.8	2.0	Thick			
ENDO SRC6048 R	60								

* L for Straight, R for Articulating, BR for Articulating with curved tip.

ENDO AFT Single Use Loading Unit

Product Code *	Staple Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Tissue Thickness	Trocar Compatibility	Number of staple rows	QTY/BOX
ENDO AFT 45TNR ENDO AFT 45TNBR	45	Tan	2/2.5/3	0.75/1/1.25	Thin	12 mm	6	6
ENDO AFT 60TNR ENDO AFT 60TNBR	60							
ENDO AFT 45PLR ENDO AFT 45PLBR	45	Purple	3/3.5/4	1.25/1.5/1.75	Regular/ Thick			
ENDO AFT 60PLR ENDO AFT 60PLBR	60							
ENDO AFT 45BKR ENDO AFT 60BKR	45 60	Black	4/4.5/5	1.75/2/2.2	Thick	15 mm		

* R for Articulating, BR for Articulating with curved tip.

iReach SRC Powered Stapler

Product Code	Shaft Length (mm)	Handle Length (mm)	Max number of reloads	QTY/BOX
ENDO SRC PSZ/PS	72	326	25	3
ENDO SRC PZ/P	162	416		
ENDO SRC PLZ/PL	262	516		



ACS-D Circular Stapler



Product Code	Color	Shaft Length	Outer Diameter (mm)	Lumen (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of Staples	Rows of Staples	QTY/BOX
ACS 21D	Gray	27 cm	21.5	12.5	3.5	1.5-1.8	16	2	3
ACS 25D	White		25.5	16	4.8	1.5-2.2	22		
ACS 28D	Blue		28.5	18.8			26		
ACS 31D	Green		31.5	22.7			30		
ACS 34D	Gold		34.5	25			32		

PPH 33D Procedure Set for Hemorrhoid and Prolapse



Product Code	Cutting Area Diameter (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
PPH 33D	24.4	3.5	0.8 ~ 1.5	32	2	3

RLS Linear Stapler



Product Code	Stapler Line (mm)	Preloaded with unit	Max Firings	QTY/BOX
RLS3035S	30	RLS3035L	8	3
RLS4535S	45	RLS4535L		
RLS6035S	60	RLS6035L		
RLS3048S	30	RLS3048L		
RLS4548S	45	RLS4548L		
RLS6048S	60	RLS6048L		

RLS Linear Stapler Cartridge

Product Code	Stapler Line (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
RLS3035L	30	Blue	3.5	1.5	11	2	6
RLS4535L	45				15		
RLS6035L	60				21		
RLS3048L	30	Green	4.8	2.0	11		
RLS4548L	45				15		
RLS6048L	60				21		

FLC Linear Cutter



Product Code	Stapler Line Length (mm)	Preloaded with cartridge	Max Firings	QTY/BOX
FLC6038S	60	FLC6038L	8	3
FLC6048S	60	FLC6048L		
FLC8038S	80	FLC8038L		
FLC8048S	80	FLC8048L		

FLC Linear Cutter Cartridge




Product Code	Stapler Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staple rows	QTY/BOX
FLC6038L	60	Blue	3.8	1.5	4	6
FLC6048L	60	Green	4.8	2.0		
FLC8038L	80	Blue	3.8	1.5		
FLC8048L	80	Green	4.8	2.0		

RLC Linear Cutter



Product Code	Stapler Line Length (mm)	Preloaded with cartridge	Max Firings	QTY/BOX
RLC6038S	60	RLC6038L	8	3
RLC6048S	60	RLC6048L		
RLC8038S	80	RLC8038L		
RLC8048S	80	RLC8048L		

RLC Linear Cutter Cartridge



Product Code	Stapler Line Length (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staple rows	QTY/BOX
RLC6038L	60	Blue	3.8	1.5	4	6
RLC6048L	60	Green	4.8	2.0		
RLC8038L	80	Blue	3.8	1.5		
RLC8048L	80	Green	4.8	2.0		

Sound Reach energy platform

Product Code	Description	QTY/BOX
CSUS8000	Generator	1
TRA6	Transducer	1
FSW2	Foot switch	1

Sound Reach Shears

Product Code	Shaft Length (cm)	QTY/ BOX
CH14PD	14	3
CH23PD	23	3
CH36PD	36	3
CH45PD	45	3



Contacts

Global Headquarter
120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25326259
F +86-22-66200915

**Medical Device Safety
Service GmbH**
Schiffgraben 41, 30175
Hanover, Germany

T +49-511-6262-8630
F +49-511-6262-8633



Importer



Distributor

