




EndoRings™

Distal Attachment



Are you getting the full picture?

- Improve your adenoma detection rate
- Compatible with all leading colonoscopes



Improved
visibility




Scope centering
during screening



Anchoring during
endoscopic therapy

The EndoRings™ device features...

- 3ayers of flexible, soft circular rings – gently flatten folds during withdrawal for improved visibility.
- Up to 50mm diameter – provides scope centering and stability of the endoscope during withdrawal and endoscopic therapy.



The EndoRings™ device easily attaches to the distal end of adult and slim colonoscopes to provide:




Improved visibility



Scope centering during screening



Anchoring during endoscopic therapy





Result from the multi-centre "CLEVER study" showed that in 126 patients, the Adenoma Detection Rate was 51% with the EndoRings versus 29% with a standard colonoscope.

The number of polyps and adenomas detected during EndoRings assisted colonoscopy was more than twice the number found during standard colonoscopy.

Data taken from "The CLEVER study" the Clinical evaluation of the EndoRings

Gentle smoothing of colonic folds allows for improved visibility!



Aquilant Endoscopy

Unit B1-B2, Pond Close, Kingsland Business Park,
Basingstoke, Hampshire, RG24 8PZ, UK
T + 44 (0) 1256 365 456 F + 44 (0) 1256 365 486
www.aquilantendoscopy.com

PART OF UDG HEALTHCARE PLC

Aquilant Endoscopy

