

Lightweight. Durable. Easy To Clean

## Enterprise Table

The Enterprise table range has been designed to meet the requirements of modern healthcare environments.


For complete flexibility they have been designed for the patient to use both in bed and when seated in a chair.



Enterprise Overbed/chair Table

**Product Code: PT10K**

**Price:**

Furniture item	
Product code	<b>PT10K</b>
Height range	710mm - 1080mm
Table top	900mm x 420mm



**Optional:**  
MBC: Plastic base cover

**Colour options available:**  
Standard; Beech  
Optional; Blue & Grey



### Features:

- 50mm twin wheel castors for easy manoeuvrability
- Gas spring assisted height adjustment
- Excellent height range of 710 – 1080mm
- Smooth, non-absorbent surfaces for ease of cleaning
- Safe working load (evenly distributed) 25kg
- Product weight 17kg