


FAVITOR

make your retraction simple

unparalleled access and visualization
increased productivity for assistants
occupies less space than metal retractors
ease of readjustment without losing focus
reduced post-op pain

RETRACTOR RINGS



C11

14.4 × 14.4 cm


Colon and Rectal Surgeries
General Surgeries
Head and Neck



C12

16.6 × 16.6 cm


Gynaecology



C13

25.4 × 14.4 cm


General Surgeries
Head and Neck



C14

18.6 × 8.3 cm

Head and Neck
Pediatrics



C15

25.3 × 11.3 cm

Gynaecology
Head and Neck



C16

32.5 × 18.3 cm

Gynaecology
Urology



C17

22 × 11.3 cm

Gynaecology
especially narrow
pelvis, obese







C18

29.1 × 18.3 cm

Gynaecology
Male Urology

ELASTIC STAYS



SKU No.	Description	Packaging
C21-8	5mm sharp hook	8/blister pack
C22-8	12mm blunt hook	8/blister pack
C23-4	7mm 2-finger sharp hook	4/blister pack
C24-4	12mm 2-finger blunt hook	4/blister pack
C31-8	5mm sharp rain drop hook	8/blister pack
C32-8	12mm blunt rain drop hook	8/blister pack
C33-4	7mm 2-finger sharp rain drop hook	4/blister pack
C34-4	12mm 2-finger blunt rain drop hook	4/blister pack