



R82 Flamingo potty frame

Comfortable potty training
and bathtime fun

R82[®] 



R82 Flamingo potty frame

We are adding a new accessory to the Flamingo family, which provides several new uses for the already highly adaptable Flamingo. Simply mount our existing Flamingo seat on the potty frame and it can be used as a freestanding device for potty training or bathtime fun. The inlay of the seat provides comfort for the user while the small size makes it fit into almost every bathing situation or room.

Features

Inlay

A good seating position is easily obtained in the comfortably shaped Flamingo seat with inlays.

Accessories

The Flamingo seat has a range of accessories.

Small frame


Simple and lightweight frame that features height adjustability

Back recline

Back recline option for the child's comfort and positioning



Versatility



Bathtime is fun and relaxing

Flamingo seat with a low potty frame solution is the perfect choice if the bathroom demands a smaller bath chair solution. The inlay and the back recline function make it very easy to be comfortable and relaxed.

Potty training

The Flamingo potty chair solution is designed to assist the child in potty training. The seat comes with an inlay and has numerous accessory options that make it easy to obtain a proper and comfortable seating position.

Frame and seat combination




Low, height adjustable frame in two different sizes



Size 1 frame fits the size 1 Flamingo seat
Size 2 frame fits the size 2 Flamingo seat
Flamingo seats available with grey or blue PUR inlay

Accessories

Numerous accessory options can be chosen for the Flamingo seat. You can explore some of the accessories below and see the complete Flamingo seat and its range of accessory options on our website.



- 1 Sides w/upholstery**
88011x
- 2 Back extension**
880118, 880119
- 3 Rod for mounting back extension**
880120
- 4 Armrests/tray**
88013x
- 5 Commode pan**
88012x
- 6 Splash guard, low**
88019x
- 7 Inlay with wider sides**
88016x
- 8 Seat inlay**
88018x



Measurements

Flamingo potty frame measurements		size 1	size 2
	Approximate total user height	990 - 1150 mm	1140 - 1300 mm
	Overall width, frame	320 mm	360 mm
	Overall Length, frame	295 mm	400 mm
	Seat width	240 mm	280 mm
	Seat depth	250 mm	290 mm
	Backrest height	350 mm	380 mm
	Seat surface height at front edge	190 - 235 mm	250 - 295 mm
	Back width	200 mm	220 mm
	Overall height	535 - 580 mm	620 - 665 mm
	Backrest angle	-90° - 30°	-90° - 30°
	Max load, total	35 kg	50 kg

Flamingo range of possibilities

Our Flamingo family is one of the most versatile products on the market with up to five configurations available. It can be used freestanding with a pan, directly over the toilet, as a shower chair and as a potty chair. The seat is also removable for direct mounting on a toilet or in a bath.



Flamingo High-low

Flamingo

Flamingo potty frame

Flamingo seat for bath tub