



Standing

R82<sup>®</sup>



User height  
700–1650 mm



Max. load  
80 kg



Environment  
Indoor only

# Gazelle PS<sup>®</sup>

Gazelle PS stands for Prone and Supine which means it can be used by children who need to be supported on their fronts or backs. It offers up to 30 degrees of hip abduction in each leg. The Gazelle PS is available in 3 sizes.

## Highlights about Gazelle PS



### Tilt vertical-horizontal

The Gazelle PS tilts from vertical to 15° past horizontal. This is an excellent position for respiratory therapy and postural drainage.



### Compatible accessories

The Gazelle accessories are compatible with all sizes, and are easily fitted and adjusted on the Gazelle frame.



### Multi-adjustable foot supports

Ball joints in foot plates allow for positioning and support for almost any ankle/foot placement.



### Manual leg abduction

Gazelle moves in an anatomically correct motion to help position the head of the femur into the acetabulum. You can adjust either one or both legs to the optimal angle.



## Standard

Multi adjustable foot support, knee support, foldable bracket for back support, removable back support, trunk safety belt, castors with lock, manual leg abduction, tilt from vertical to horizontal, adjustable trunk or chest support.

## Options

### Size

- Size 1
- Size 2
- Size 3

### Frame colour



Yellow for size 1



Green for size 2



Blue for size 3



As standard the Gazelle PS is supplied with chest and foot straps

## Measurements

		1	2	3
	User height: floor to axil	450-950 mm	600-1050 mm	700-1250 mm
	Max. load	40 kg	60 kg	80 kg
	Angle, center bar	90 - -15°	90 - -15°	90 - -15°
	Length	760 mm	880 mm	990 mm
	Overall width	480 mm	580 mm	650 mm
	Overall height	760 mm	800 mm	900 mm
	Total mass	20 kg	22 kg	24 kg

Further measurements at our website

## Accessories



**Head support system**  
88124-xx



**Mounting for head support**  
88130



**Facial aperture**  
88117-x



**Anat. head support**  
88119-x



**Side supports**  
88055-x



**Side supports, interior**  
88075-x



**Side supports, exterior**  
88065-x



**Trunk support**  
88106-x



**Angling of lower leg**  
88165-x



**Strap for side supports**  
88250-x



**Tray with bowl**  
88113-x



**Tray**  
88108-x

**Padding for tray**  
862135-x

For more accessories - visit our website