



iReach Magnum

Powered Endoscopic
Linear Stapler

55°

55° articulation



Powered
12 Firings



Superior
Stability



Superior
Compression



New blade in
every cartridge

Product Code	Staple Line Length	Cartridge Color	Open Staple Height	Closed Staple Height	Tissue Thickness	Number of staple rows	QTY/ BOX
REC 45 GRA	45mm	Gray	2.0 mm	0.75 mm	Extra Thin	6	12
REC 60 GRA	60mm						
REC 45 WHT	45mm	White	2.5 mm	1.0 mm	Vascular		
REC 60 WHT	60mm						
REC 45 BLU	45mm	Blue	3.5 mm	1.5 mm	Thin		
REC 60 BLU	60mm						
REC 45 GLD	45mm	Gold	3.8 mm	1.75 mm	Medium		
REC 60 GLD	60mm						
REC 45 GRN	45mm	Green	4.1 mm	2.0 mm	Thick		
REC 60 GRN	60mm						
REC 45 BLK	45mm	Black	4.4 mm	2.2 mm	Extra Thick		
REC 60 BLK	60mm						

Contacts

Global Headquarter
120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25326259
F +86-22-66200915

**Medical Device Safety
Service GmbH**
Schiffgraben 41, 30175
Hanover, Germany

T +49-511-6262-8630
F +49-511-6262-8633



Importer



Distributor

