



ORDERING INFORMATION

Fax order to Customer Services on 01 885 5558

	BSN Code	Option	Size	Length	Colour	PIP Code	NHSSC Code*
	CCL1	Knitted Band	Small	Regular	Beige	382-1196	EGD7879
		Knitted Band	Small	Long	Beige	325-7011	EGD340
		Knitted Band	Medium	Regular	Beige	382-1188	EGD7880
		Knitted Band	Medium	Long	Beige	325-7037	EGD341
		Knitted Band	Large	Regular	Beige	382-1170	EGD7881
		Knitted Band	Large	Long	Beige	325-7029	EGD342
		Silicone Band	Small	Regular	Beige	382-1238	EGD7882
		Silicone Band	Small	Long	Beige	325-6955	EGD346
		Silicone Band	Medium	Regular	Beige	382-1220	EGD7883
		Silicone Band	Medium	Long	Beige	325-6971	EGD347
	CCL2	Knitted Band	Small	Regular	Beige	382-1261	EGD7885
		Knitted Band	Small	Long	Beige	325-6989	EGD343
		Knitted Band	Medium	Regular	Beige	382-1204	EGD7886
		Knitted Band	Medium	Long	Beige	325-7003	EGD344
		Knitted Band	Large	Regular	Beige	382-1246	EGD7887
		Knitted Band	Large	Long	Beige	325-6997	EGD345
		Silicone Band	Small	Regular	Beige	382-1287	EGD7888
		Silicone Band	Small	Long	Beige	325-6898	EGD7891
		Silicone Band	Medium	Regular	Beige	382-1279	EGD7889
		Silicone Band	Medium	Long	Beige	325-6914	EGD350
		Silicone Band	Large	Regular	Beige	382-1253	EGD7890
		Silicone Band	Large	Long	Beige	325-6906	EGD351
		Knitted Band	Small	Long	Black	363-7568	
		Knitted Band	Medium	Long	Black	363-7576	
		Knitted Band	Large	Long	Black	363-7550	

GAUNTLET ORDERING INFORMATION

	BSN Code	Option	Size	Length	Colour	PIP Code
	CCL1	With Thumb	Small	-	Beige	325-6823
		With Thumb	Medium	-	Beige	325-6831
		With Thumb	Large	-	Beige	325-6922
	CCL2	With Thumb	Small	-	Beige	325-6773
		With Thumb	Medium	-	Beige	325-6781
		With Thumb	Large	-	Beige	325-6799
		With Thumb	Small	-	Black	363-7543
		With Thumb	Medium	-	Black	363-7535
		With Thumb	Large	-	Black	363-7527

PLEASE TURN OVER FOR COMBINED ARMSLEEVE ORDERING CODES AND MEASURING GUIDES

JOBST® BELLA™ LITE



READY-TO-WEAR UPPER LIMB COMPRESSION GARMENTS COMBINED

ARMSLEEVE ORDERING INFORMATION

Fax order to Customer Services on 01 885 5558

HC21
Healthcare

JOBST
Comfort, Health and Style!


	BSN Code	Option	Size	Length	Colour	PIP Code
	77678-00	Knitted Band	Small	Regular	Beige	398-9571
	77678-03	Knitted Band	Small	Long	Beige	398-9605
	77678-01	Knitted Band	Medium	Regular	Beige	398-9589
	77678-04	Knitted Band	Medium	Long	Beige	398-9613
	77678-02	Knitted Band	Large	Regular	Beige	398-9597
	77678-05	Knitted Band	Large	Long	Beige	398-9621
	77678-06	Silicone Band	Small	Regular	Beige	398-9639
	77678-09	Silicone Band	Small	Long	Beige	398-9662
	77678-07	Silicone Band	Medium	Regular	Beige	398-9647
	77678-10	Silicone Band	Medium	Long	Beige	398-9670
	77678-08	Silicone Band	Large	Regular	Beige	398-9654
77678-11	Silicone Band	Large	Long	Beige	398-9688	
	77680-00	Knitted Band	Small	Regular	Beige	398-9696
	77680-03	Knitted Band	Small	Long	Beige	398-9720
	77680-01	Knitted Band	Medium	Regular	Beige	398-9704
	77680-04	Knitted Band	Medium	Long	Beige	398-9738
	77680-02	Knitted Band	Large	Regular	Beige	398-9712
	77680-05	Knitted Band	Large	Long	Beige	398-9746
	77680-06	Silicone Band	Small	Regular	Beige	398-9753
	77680-09	Silicone Band	Small	Long	Beige	398-9787
	77680-07	Silicone Band	Medium	Regular	Beige	398-9761
	77680-10	Silicone Band	Medium	Long	Beige	398-9795
	77680-08	Silicone Band	Large	Regular	Beige	398-9779
77680-11	Silicone Band	Large	Long	Beige	398-9803	

MEASURING GUIDE

	Armsleeve			Gauntlet	
	Wrist	Mid-Lower	Mid-Upper	Palm	Wrist
Small	14cm - 20cm	18cm - 23cm	22cm - 40cm	15cm - 19cm	14cm - 20cm
Medium	16cm - 22cm	23cm - 29cm	27cm - 45cm	19cm - 22cm	16cm - 22cm
Large	18cm - 24cm	29cm - 34cm	33cm - 51cm	22cm - 26cm	18cm - 24cm
	Combined Armsleeve			Length	
	Palm	Mid-Lower	Mid-Upper	Length	
Small	15cm-19cm	18cm-23cm	22cm-40cm	Regular up to (≤) 43 cm	
Medium	19cm-22cm	23cm-29cm	27cm-45cm	Long over (>) 43 cm	
Large	22cm-26cm	29cm-34cm	33cm-51cm	Note: Length measurement taken on inside of arm	

PLEASE TURN OVER FOR ARMSLEEVE AND GAUNTLET ORDERING CODES.

Colours Actual garment colour may vary from printed swatch sample



Customer Service: 01 885 5551 Compression Therapy Helpline: +44 800 389 8424
Email: info@hc21.ie www.healthcare21.eu

