



Comfort, Health and Style!

JOBST® Bella™ Lite

Comfortable and effective

Ready-to-wear,
circular-knit, upper limb
compression garments



THERAPIES. HAND IN HAND.



TWO LENGTHS FOR OPTIMAL FIT

JOBST BELLA LITE IS DESIGNED TO TREAT THE SYMPTOMS OF MILD TO MODERATE LYMPHOEDEMA.

These comfortable garments are available as an armsleeve or a gauntlet that can be worn together or separately and a combined armsleeve.

Two compression classes

- Compression class 1 (15-20mmHg) and 2 (20-30mmHg)
- Effective graduated compression suitable for treating mild to moderate levels of lymphoedema

Style choice

- Three styles, armsleeve, gauntlet with thumb and combined armsleeve and gauntlet
- Armsleeve available with knitted or dotted silicone band allowing for individual patient choice
- Two colours, beige or black*

Three sizes and two lengths

- Fits a wide range of arm and hand measurements
- Only 2-3 measurements required per fitting
- Two armsleeve lengths for optimal fit

Soft fine-knit construction

- Garments are soft, durable and air permeable ensuring they are comfortable to wear and easy to put on
- Smooth, flat bands offer reduced compression at the wrist to provide an effective and comfortable fit when worn together
- Latex-free formulation ensures reduced risk of patient sensitivity or skin reaction

Available on prescription

- All garments are fully reimbursable

For ordering details, please contact your **BSN medical** Account Manager or visit www.bsnmedical.co.uk.



*CCL 2 knitted band armsleeve in long and gauntlet only