

# JOBST® CONFIDENCE

## Custom Fit Order Form



TO ORDER PLEASE EMAIL: info@hc21.group

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_

Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_

Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_



<input type="checkbox"/> <b>JOBST® Confidence</b>		
<b>Compression Class (RAL)</b>	<b>Quantity</b>	
<input type="checkbox"/> CCL 1 (18-21mmHg)	Left	<input type="text"/>
<input type="checkbox"/> CCL 2 (23-32mmHg)	Right	<input type="text"/>
<input type="checkbox"/> CCL 3 (34-46mmHg)		
<b>Style and Options</b>		
<input type="checkbox"/> AD Knee high		
<input type="checkbox"/> Slipform	<input type="checkbox"/> Ankle functional zone	<input type="checkbox"/> Hallux valgus (bunion relief zone)*
<input type="checkbox"/> Decor-line	<input type="checkbox"/> Initials (max 2 letters) <input type="text"/>	
<input type="checkbox"/> AG Thigh high		
<input checked="" type="checkbox"/> Slipform	<input type="checkbox"/> Ankle functional zone	<input type="checkbox"/> Knee functional zone
<input type="checkbox"/> Decor-line	<input type="checkbox"/> Initials (max 2 letters) <input type="text"/>	<input type="checkbox"/> Hallux valgus (bunion relief zone)*
<b>Colour</b>	<b>Top band options</b>	<b>Knee high A-D</b>
<input type="checkbox"/> Beige		<input type="checkbox"/> Black
<input type="checkbox"/> Caramel		<input type="checkbox"/> Red heather
<input type="checkbox"/> Anthracite heather	<input type="checkbox"/> Jeans heather	<b>Thigh high A-G</b>
Note: the functional zone for red, anthracite and jeans heather coloured garments is more visible		<input checked="" type="checkbox"/> 5cm silicone top band
<b>Remarks</b>		

### Left Leg Length Measurements (cms)

a-K \_\_\_\_\_

**Required for thigh high garment**

a-G \_\_\_\_\_ G \_\_\_\_\_

a-F \_\_\_\_\_ F \_\_\_\_\_

a-E \_\_\_\_\_ E \_\_\_\_\_


**a-D** \_\_\_\_\_ **D** \_\_\_\_\_  
(Knee high)

a-C \_\_\_\_\_ C \_\_\_\_\_

a-B' \_\_\_\_\_ B' \_\_\_\_\_

a-B \_\_\_\_\_ B \_\_\_\_\_

### Circumference Measurements (cms)



\*a = floor / base of measuring board

### Right Leg Length Measurements (cms)

a-K \_\_\_\_\_

**Required for thigh high garment**

G \_\_\_\_\_ a-G \_\_\_\_\_

F \_\_\_\_\_ a-F \_\_\_\_\_

E \_\_\_\_\_ a-E \_\_\_\_\_


**D** \_\_\_\_\_ **a-D** \_\_\_\_\_  
(Knee high)

C \_\_\_\_\_ a-C \_\_\_\_\_

B' \_\_\_\_\_ a-B' \_\_\_\_\_

B \_\_\_\_\_ a-B \_\_\_\_\_

### Left Foot Slant



### Foot Length - Left

For open toe \_\_\_\_\_ cm


For closed toe \_\_\_\_\_ cm  
(longest toe)

### Foot Length - Right

For open toe \_\_\_\_\_ cm

For closed toe \_\_\_\_\_ cm  
(longest toe)

### Right Foot Slant



Note: medial = inside  
lateral = outside

\*Only available with slant foot

# JOBST® CONFIDENCE

## Ordering Information

Compression Class / Option    Code    Pieces

### JOBST® Confidence Knee High



#### Basic style

CCL 1 (18-21mmHg)	L1-16-04	1
CCL 2 (23-32mmHg)	L2-17-04	1
CCL 3 (34-46mmHg)	L3-18-04	1

#### Options

Closed toe	L-A001C	1
SoftFit	L-A016C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1

Compression Class / Option    Code    Pieces

### JOBST® Confidence Thigh High



#### Basic style

CCL 1 (18-21mmHg)	L1-16-06	1
CCL 2 (23-32mmHg)	L2-17-06	1
CCL 3 (34-46mmHg)	L3-18-06	1

#### Options

Closed toe	L-A001C	1
5cm silicone band	L-A004C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1
Knee functional zone	L-A015C	1



## Options



#### Ankle Functional Zone

Circumferential functional zone allows greater movement and flexibility at the ankle



#### Knee Functional Zone

Circumferential functional zone allows greater movement and flexibility at the knee



#### Slipform

Slipform provides a higher shape to the outside of the garment, at the top, conforming to the limb shape



#### Deco Line

A decorative seam line can be added to the front of the garment on request



#### Hallux Valgus (bunion relief zone)\*

Creates a looser knitting pattern around the ball of the foot to provide relief to painful bunions



#### SoftFit Top Band

Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit



#### Initials

Patient initials can be knitted into the garment below the JOBST logo for a personal touch