


JOBST® Farrow® Glove

Date: _____ Purchase Order No.: _____ Patient Name: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____

Ready-to-wear order form and prescription request

For quotes please email: quotes@hc21.group

To order please email Customer Services on: info@hc21.group



Measuring Instructions

1. Measure the **A** circumference of the palm at the base of the fingers
2. Measure the **B** circumference around the wrist
3. Determine the size (**1-8**)

JOBST® Farrow® Glove

Size: 1 2 3 4 5 6 7 8

Colour: Tan Black

Product Code: -

Quantity Required:

Repeat prescription required every months

Size	1	2	3	4	5	6	7	8
A	16.6 - 18.0cm	16.6 - 18.0cm	18.0 - 19.5cm	18.0 - 19.5cm	18.0 - 19.5cm	19.5 - 21.5cm	19.5 - 21.5cm	21.5 - 23.1cm
B	15.6 - 16.7cm	16.7 - 18.1cm	15.9 - 17.0cm	17.0 - 18.5cm	18.5 - 20.4cm	17.2 - 18.7cm	18.7 - 20.6cm	18.3 - 20.1cm

Order Codes

Colour	Compression Level	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8
		Product Code	Product Code	Product Code	Product Code	Product Code	Product Code	Product Code	Product Code
Tan	15-20 mmHg	75463-00	75463-01	75463-02	75463-03	75463-04	75463-05	75463-06	75463-07
	20-30 mmHg	75924-00	75924-01	75924-02	75924-03	75924-04	75924-05	75924-06	75924-07
Black	15-20 mmHg	75463-10	75463-11	75463-12	75463-13	75463-14	75463-15	75463-16	75463-17
	20-30 mmHg	75924-10	75924-11	75924-12	75924-13	75924-14	75924-15	75924-16	75924-17

Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for codefinder.jobst.com on your mobile device.



By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, and the EU General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit [Group Data Protection Statement - Healthcare21](#)