


# JOBST® FarrowWrap® 4000 A-D Below Knee

Ready-to-wear order form and prescription request

For quotes please email: [quotes@hc21.group](mailto:quotes@hc21.group)  
 To order please email Customer Services on: [info@hc21.group](mailto:info@hc21.group)

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_

## Legpiece



Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

## Measuring instructions

### Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **B** (narrowest part of the ankle) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-L**) & length (**Regular or Tall**) using chart below

When ordering JOBST FarrowWrap 4000, 1 x JOBST FarrowHybrid (Liners) wide is also provided. If additional or alternative JOBST FarrowHybrid (Liners) sizes are required, please order separately

### JOBST FarrowHybrid (Liners)

1. Measure **A1** circumference
2. Record **shoe size**
3. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

Do not apply JOBST FarrowHybrid (Liners) if calf circumference exceeds 60cm. When a large JOBST FarrowWrap 4000 is ordered, an additional non-compressive liner will be provided to accommodate a calf circumference greater than 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance

### Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



**JOBST® FarrowWrap® 4000 Legpiece**

**Size:**  X Small  Small  Medium  Large

**Length**  **B-D**  Regular (30-33 cm)  Tall (34-37 cm)

**Colour:**  Tan  Black

Product Code:  -

Quantity Required:

Repeat prescription required every  months

**JOBST® FarrowHybrid® (Liners)**

**Size:**  Small (2-5)  Medium (6-8)  Large (9-11)  X Large (12-13) (shoe size)

**Width**  **A1**  Regular (20-27cm)  Wide (28-35cm)

**Colour:**  Tan  Black (available in Wide variant only)

Product Code:  -

Quantity Required:

Repeat prescription required every  months

	X Small	Small	Medium	Large
<b>JOBST FarrowWrap 4000 Legpiece</b>				
<b>C</b> Calf	25-38cm	32-48cm	37-56cm	43-67cm
<b>B</b> Ankle	17-23cm	19-28cm	25-37cm	29-43cm
<b>B-D</b> Regular	30-33cm	30-33cm	30-33cm	30-33cm
<b>B-D</b> Tall	34-37cm	34-37cm	34-37cm	34-37cm

NOTE: for **B-D** length measure posteriorly (back of leg) and follow contours

# JOBST FarrowWrap 4000 Order Codes

## JOBST® FarrowWrap® 4000 Legpiece (30-40mmHg)

For patients with **mild to moderate** oedema, with or without a leg ulcer

Product Code	Size	Length	Colour
76661-00	Extra small	Regular	Beige
76661-08	Extra small	Regular	Black
76661-04	Extra small	Tall	Beige
76661-12	Extra small	Tall	Black
76661-01	Small	Regular	Beige
76661-09	Small	Regular	Black
76661-05	Small	Tall	Beige
76661-13	Small	Tall	Black
76661-02	Medium	Regular	Beige
76661-10	Medium	Regular	Black
76661-06	Medium	Tall	Beige
76661-14	Medium	Tall	Black
76661-03	Large	Regular	Beige
76661-11	Large	Regular	Black
76661-07	Large	Tall	Beige
76661-15	Large	Tall	Black



## JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour
76664-00	Small - Regular (shoe size 2-5)	Tan
76664-01	Medium - Regular (shoe size 6-8)	Tan
76664-02	Large - Regular (shoe size 9-11)	Tan
76664-03	Extra Large - Regular (shoe size 12-13)	Tan
76664-12	Small - Wide (shoe size 2-5)	Tan
76664-13	Medium - Wide (shoe size 6-8)	Tan
76664-14	Large - Wide (shoe size 9-11)	Tan
76664-15	Extra Large - Wide (shoe size 12-13)	Tan
76664-16	Small - Wide (shoe size 2-5)	Black
76664-17	Medium - Wide (shoe size 6-8)	Black
76664-18	Large - Wide (shoe size 9-11)	Black
76664-19	Extra Large - Wide (shoe size 12-13)	Black



Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance

### Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for **codefinder**. [jobst.com](http://jobst.com) on your mobile device.



Scan to access HC21 Compression Products available on PCRS



SCAN ME

**HC21**  
Healthcare

info@hc21.group  
00 353 214 875 055  
healthcare21.eu

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018 and the (EU) General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit [www.healthcare21.eu](http://www.healthcare21.eu)