

# JOBST® FarrowWrap®

## A-D Below Knee / A-G Thigh High

Ready-to-Wear and Custom Fit Order Form.




TO ORDER PLEASE EMAIL: [info@hc21.group](mailto:info@hc21.group)

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Email: \_\_\_\_\_  
 \_\_\_\_\_ Invoice Address: \_\_\_\_\_  
 \_\_\_\_\_



Please measure all points for a JOBST® FarrowWrap® custom-fit product



**E-G**

Left: \_\_\_\_\_

Right: \_\_\_\_\_

**D-E**

Left: \_\_\_\_\_

Right: \_\_\_\_\_

Required for thighpiece

**A-D**

Left: \_\_\_\_\_

Right: \_\_\_\_\_

**Note:**  
For A-D and D-E, measure posteriorly and follow leg contour.

**Note:**  
For E-G measure medially and follow contours.

Custom-fit

Brand Compression	Quantity:	LEFT	RIGHT
Lite	<input type="text"/>	<input type="text"/>	<input type="text"/>
Classic	<input type="text"/>	<input type="text"/>	<input type="text"/>
Strong	<input type="text"/>	<input type="text"/>	<input type="text"/>

**Remarks**

Please measure all points below for a JOBST® FarrowWrap® ready-to-wear product

	X Small	Small	Medium	Large	X Large	XX Large
<b>Thighpiece (including kneepiece)</b>						
<b>G</b> Groin	60-70 cm	65-75 cm	70-80 cm	75-90 cm	85-100 cm	-
<b>E</b> Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	-
<b>E-G</b> Short	16-20 cm	16-20 cm	16-20 cm	16-20 cm	16-20 cm	-
<b>E-G</b> Regular	21-26 cm	21-26 cm	21-26 cm	21-26 cm	21-26 cm	-
<b>E-G</b> Tall	27-31 cm	27-31 cm	27-31 cm	27-31 cm	27-31 cm	-
NOTE: for <b>E-G</b> measure medially and follow contours						
<b>Kneepiece (seperate)</b>						
<b>E</b> Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	81-87 cm
<b>C</b> Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	63-73 cm
<b>Legpiece</b>						
<b>C</b> Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	-
<b>B</b> Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm	-
<b>A-D</b> Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm	-
<b>A-D</b> Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm	-
NOTE: for <b>A-D</b> measure posteriorly and follow contours						
<b>Footpiece</b>						
<b>A¹</b> Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm	-
<b>X</b> Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm	-
<b>X</b> Tall	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm	-
NOTE: for <b>X</b> measure from the heel to the base of the little toe						

**Remarks**

Note: 3cm allowance on either side of measurement.

Please use ready-to-wear product code required found on the back of this form

**Ready-to-wear**

	Quantity Left	Product Code	Quantity Right	Product Code
<b>Thigh Piece</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Leg Piece</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Foot Piece</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

A pair of appropriately sized liners will be provided free of charge with JOBST® FarrowWrap® lower limb.

# JOBST® FarrowWrap® ready-to-wear direct order codes

## JOBST® FarrowWrap® Lite (20-30mmHg)

Code	Size	Length
<b>Thighpiece (including kneepiece)</b>		
76655-05	Extra Small	Short
76655-00	Extra Small	Regular
76655-10	Extra Small	Tall
76655-06	Small	Short
76655-01	Small	Regular
76655-11	Small	Tall
76655-07	Medium	Short
76655-02	Medium	Regular
76655-12	Medium	Tall
76655-08	Large	Short
76655-03	Large	Regular
76655-13	Large	Tall
76655-09	Extra Large	Short
76655-04	Extra Large	Regular
76655-14	Extra Large	Tall
<b>Kneepiece (seperate)*</b>		
73367-00	Extra Small	-
73367-01	Small	-
73367-02	Medium	-
73367-03	Large	-
73367-04	Extra Large	-
<b>Legpiece</b>		
76654-00	Extra Small	Regular
76654-05	Extra Small	Tall
76654-01	Small	Regular
76654-06	Small	Tall
76654-02	Medium	Regular
76654-07	Medium	Tall
76654-03	Large	Regular
76654-08	Large	Tall
76654-04	Extra Large	Regular
76654-09	Extra Large	Tall
<b>Footpiece</b>		
76652-05	Extra Small	Regular
76652-00	Extra Small	Long
76652-06	Small	Regular
76652-01	Small	Long
76652-07	Medium	Regular
76652-02	Medium	Long
76652-08	Large	Regular
76652-03	Large	Long
76652-09	Extra Large	Regular
76652-04	Extra Large	Long



## JOBST® FarrowWrap® Strong (30-40mmHg)

Code	Size	Length
<b>Thighpiece (including kneepiece)</b>		
76649-05	Extra Small	Short
76649-00	Extra Small	Regular
76649-10	Extra Small	Tall
76649-06	Small	Short
76649-01	Small	Regular
76649-11	Small	Tall
76649-07	Medium	Short
76649-02	Medium	Regular
76649-12	Medium	Tall
76649-08	Large	Short
76649-03	Large	Regular
76649-13	Large	Tall
76649-09	Extra Large	Short
76649-04	Extra Large	Regular
76649-14	Extra Large	Tall
<b>Kneepiece (seperate)*</b>		
73368-00	Extra Small	-
73368-01	Small	-
73368-02	Medium	-
73368-03	Large	-
73368-04	Extra Large	-
73368-05	Extra Extra Large	-
<b>Legpiece</b>		
76648-00	Extra Small	Regular
76648-05	Extra Small	Tall
76648-01	Small	Regular
76648-06	Small	Tall
76648-02	Medium	Regular
76648-07	Medium	Tall
76648-03	Large	Regular
76648-08	Large	Tall
76648-04	Extra Large	Regular
76648-09	Extra Large	Tall
<b>Footpiece</b>		
76646-04	Extra Small	Regular
76646-09	Extra Small	Long
76646-05	Small	Regular
76646-00	Small	Long
76646-06	Medium	Regular
76646-01	Medium	Long
76646-07	Large	Regular
76646-02	Large	Long
76646-08	Extra Large	Regular
76646-03	Extra Large	Long



## JOBST® FarrowWrap® Classic (30-40mmHg)

Code	Size	Length
<b>Thighpiece (including kneepiece)</b>		
76658-06	Extra Small	Short
76658-00	Extra Small	Regular
76658-12	Extra Small	Tall
76658-07	Small	Short
76658-01	Small	Regular
76658-13	Small	Tall
76658-08	Medium	Short
76658-02	Medium	Regular
76658-14	Medium	Tall
76658-09	Large	Short
76658-03	Large	Regular
76658-15	Large	Tall
76658-10	Extra Large	Short
76658-04	Extra Large	Regular
76658-16	Extra Large	Tall
<b>Legpiece</b>		
76657-00	Extra Small	Regular
76657-05	Extra Small	Tall
76657-01	Small	Regular
76657-06	Small	Tall
76657-02	Medium	Regular
76657-07	Medium	Tall
76657-03	Large	Regular
76657-08	Large	Tall
76657-04	Extra Large	Regular
76657-09	Extra Large	Tall
<b>Footpiece</b>		
76656-05	Extra Small	Regular
76656-00	Extra Small	Long
76656-06	Small	Regular
76656-01	Small	Long
76656-07	Medium	Regular
76656-02	Medium	Long
76656-08	Large	Regular
76656-03	Large	Long
76656-09	Extra Large	Regular
76656-04	Extra Large	Long

### Garment Grip Roll\*

Code	Size
<b>Garment Grip Roll</b>	
76672-00	14cm x 80cm
76672-01	14cm x 120cm
76672-02	14cm x 550cm
76672-03	20cm x 120cm
76672-04	30cm x 150cm
76672-05	10cm x 60cm

## JOBST® FarrowHybrid® - supplying compression to the foot and ankle

Code	Size	Colour
<b>ADI (20-30mmHg)</b>		
76664-00	Small (shoe size 2-5)	Beige
76664-01	Medium (shoe size 6-8)	Beige
76664-02	Large (shoe size 9-11)	Beige
76664-03	Extra large (shoe size 12-13)	Beige
76664-12	Small - wide (shoe size 2-5)	Beige
76664-13	Medium - wide (shoe size 6-8)	Beige
76664-14	Large - wide (shoe size 9-11)	Beige
76664-15	Extra large - wide (shoe size 12-13)	Beige
76664-16	Small - wide (shoe size 2-5)	Black
76664-17	Medium - wide (shoe size 6-8)	Black
76664-18	Large - wide (shoe size 9-11)	Black
76664-19	Extra large - wide (shoe size 12-13)	Black



\*Do not apply JOBST® FarrowHybrid if calf circumference exceeds 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance.

Wide = A1 Measurement > 28cm in circumference

### A-D Below Knee Liners\*

Code	Size
<b>Beige</b>	
76666-21	Small (15-58cm)
76666-22	Large (52-70cm)
<b>TG Soft</b>	
76669-00	Small (<40cm)
76669-02	Medium (40-70cm)
76669-04	Large (70-120cm)

### A-G Thigh High Liners\*

Code	Size
<b>TG Soft</b>	
76670-02	Medium (40-70cm)
76670-04	Large (70-120cm)

### Foam Squares Pack\*

Code	Size
<b>Foam Squares Pack</b>	
76675-01	Small (8x8cm)
76675-02	Large (10x10cm)

NB The thighpiece comes with kneepiece included, however, the kneepiece can be ordered separately.

