

# JOBST® FarrowWrap® A-G Thigh High

Ready-to-wear order form and prescription request


For quotes please email: [quotes@hc21.group](mailto:quotes@hc21.group)

To order please email Customer Services on: [info@hc21.group](mailto:info@hc21.group)

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_

Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_

Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_



Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

## Measuring instructions

### Thighpiece

1. Measure **E** (5cm above the patella) and **G** (2cm below gluteal fold) circumferences
2. Measure length medially (inside leg) following the contour from **E** to **G**
3. Determine size (**XS-XL**) & length (**Short, Regular or Tall**) using chart below  
A knee piece & a pair of appropriately sized TG Soft Liners will be provided free of charge with each JOBST FarrowWrap thighpiece

### Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **A** (floor) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-XL**) & length (**Regular or Tall**) using chart below.  
A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap legpiece

### Footpiece

1. Measure **A1** circumference
2. Measure length **X** from heel to base of little toe
3. Determine size (**XS-XL**) & length (**Regular or Long**) using chart below

### JOBST FarrowHybrid (Liners)

1. JOBST FarrowHybrid (Liners) can be used instead of a footpiece, providing calf circumference does not exceed 60cm
2. Measure **A1** circumference
3. Record shoe size
4. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

### Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



	X Small	Small	Medium	Large	X Large
<b>Thighpiece (including kneepiece)</b>					
<b>G</b> Groin	60-70cm	65-75cm	70-80cm	75-90cm	85-100cm
<b>E</b> Above knee	47-52cm	53-59cm	60-66cm	67-73cm	74-80cm
<b>E-G</b> Short	16-20cm	16-20cm	16-20cm	16-20cm	16-20cm
<b>E-G</b> Regular	21-26cm	21-26cm	21-26cm	21-26cm	21-26cm
<b>E-G</b> Tall	27-31cm	27-31cm	27-31cm	27-31cm	27-31cm
NOTE: for <b>E-G</b> length measure medially (inside leg) and follow contours					
<b>Legpiece</b>					
<b>C</b> Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
<b>B</b> Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
<b>A-D</b> Regular	33-37cm	35-39cm	37-41cm	39-43cm	39-43cm
<b>A-D</b> Tall	38-41cm	40-43cm	42-45cm	44-47cm	44-47cm
NOTE: for <b>A-D</b> length measure posteriorly (back of leg) and follow contours					
<b>Footpiece</b>					
<b>A1</b> Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
<b>X</b> Regular	16-17cm	17-18cm	19-20cm	20-21cm	22-23cm
<b>X</b> Long	18-19cm	19-20cm	21-22cm	22-23cm	24-25cm
NOTE: for <b>X</b> length measure from the heel to the base of the little toe					

Note: 3cm allowance on either side of measurement

### Thighpiece

**Size:**  X Small  Small  Medium  Large  X Large

**Length **E-G****  Short  Regular  Tall

**Fabric:**  Lite  Strong  Classic

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

### Legpiece

**Size:**  X Small  Small  Medium  Large  X Large

**Length **A-D****  Regular  Tall

**Fabric:**  Lite  Strong  Classic

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

### Footpiece

**Size:**  X Small  Small  Medium  Large  X Large

**Length **X****  Regular  Long

**Fabric:**  Lite  Strong  Classic

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

### JOBST® FarrowHybrid® (Liners)

**Size:**  Small (2-5)  Medium (6-8)  Large (9-11)  X Large (12-13) (shoe size)

**Width **A1****  Regular (20-27cm)  Wide (28-35cm)

**Colour:**  Tan  Black (available in Wide variant only)

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

# JOBST FarrowWrap Order Codes

## JOBST® FarrowWrap® Lite (20-30mmHg)

For patients with **mild to moderate** oedema caused by venous or lymphatic conditions, where reduced compression is required

Product Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76655-05	Extra Small	Short	Beige	360-6084
76655-00	Extra Small	Regular	Beige	360-6134
76655-10	Extra Small	Tall	Beige	360-6183
76655-06	Small	Short	Beige	360-6092
76655-01	Small	Regular	Beige	360-6142
76655-11	Small	Tall	Beige	360-6191
76655-07	Medium	Short	Beige	360-6100
76655-02	Medium	Regular	Beige	360-6159
76655-12	Medium	Tall	Beige	360-6209
76655-08	Large	Short	Beige	360-6118
76655-03	Large	Regular	Beige	360-6167
76655-13	Large	Tall	Beige	360-6217
76655-09	Extra Large	Short	Beige	360-6126
76655-04	Extra Large	Regular	Beige	360-6175
76655-14	Extra Large	Tall	Beige	360-6225
<b>Legpiece</b>				
76654-00	Extra Small	Regular	Beige	360-6233
76654-05	Extra Small	Tall	Beige	360-6282
76654-01	Small	Regular	Beige	360-6241
76654-06	Small	Tall	Beige	360-6290
76654-02	Medium	Regular	Beige	360-6258
76654-07	Medium	Tall	Beige	360-6308
76654-03	Large	Regular	Beige	360-6266
76654-08	Large	Tall	Beige	360-6316
76654-04	Extra Large	Regular	Beige	360-6274
76654-09	Extra Large	Tall	Beige	360-6324
<b>Footpiece</b>				
76652-05	Extra Small	Regular	Beige	360-6332
76652-00	Extra Small	Long	Beige	360-6381
76652-06	Small	Regular	Beige	360-6340
76652-01	Small	Long	Beige	360-6399
76652-07	Medium	Regular	Beige	360-6357
76652-02	Medium	Long	Beige	360-6407
76652-08	Large	Regular	Beige	360-6365
76652-03	Large	Long	Beige	360-6415
76652-09	Extra Large	Regular	Beige	360-6373
76652-04	Extra Large	Long	Beige	360-6423



## JOBST® FarrowWrap® Strong (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for skin folds and shape distortion

Product Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76649-05	Extra Small	Short	Beige	382-0222
76649-00	Extra Small	Regular	Beige	382-0172
76649-10	Extra Small	Tall	Beige	382-0123
76649-06	Small	Short	Beige	382-0255
76649-01	Small	Regular	Beige	382-0206
76649-11	Small	Tall	Beige	382-0156
76649-07	Medium	Short	Beige	382-0248
76649-02	Medium	Regular	Beige	382-0198
76649-12	Medium	Tall	Beige	382-0149
76649-08	Large	Short	Beige	382-0230
76649-03	Large	Regular	Beige	382-0180
76649-13	Large	Tall	Beige	382-0131
76649-09	Extra Large	Short	Beige	382-0214
76649-04	Extra Large	Regular	Beige	382-0164
76649-14	Extra Large	Tall	Beige	382-0115
<b>Legpiece</b>				
76648-00	Extra Small	Regular	Beige	382-0073
76648-05	Extra Small	Tall	Beige	382-0024
76648-01	Small	Regular	Beige	382-0107
76648-06	Small	Tall	Beige	382-0057
76648-02	Medium	Regular	Beige	382-0099
76648-07	Medium	Tall	Beige	382-0040
76648-03	Large	Regular	Beige	382-0081
76648-08	Large	Tall	Beige	382-0032
76648-04	Extra Large	Regular	Beige	382-0065
76648-09	Extra Large	Tall	Beige	382-0016
<b>Footpiece</b>				
76646-04	Extra Small	Regular	Beige	381-9976
76646-09	Extra Small	Long	Beige	381-9927
76646-05	Small	Regular	Beige	382-0008
76646-00	Small	Long	Beige	381-9950
76646-06	Medium	Regular	Beige	381-9992
76646-01	Medium	Long	Beige	381-9943
76646-07	Large	Regular	Beige	381-9984
76646-02	Large	Long	Beige	381-9935
76646-08	Extra Large	Regular	Beige	381-9968
76646-03	Extra Large	Long	Beige	381-9919



## JOBST® FarrowWrap® Classic (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for irregular shaped limbs and deep skin folds

Product Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76658-06	Extra Small	Short	Beige	360-6431
76658-00	Extra Small	Regular	Beige	360-6480
76658-12	Extra Small	Tall	Beige	360-6530
76658-07	Small	Short	Beige	360-6449
76658-01	Small	Regular	Beige	360-6498
76658-13	Small	Tall	Beige	360-6548
76658-08	Medium	Short	Beige	360-6456
76658-02	Medium	Regular	Beige	360-6506
76658-14	Medium	Tall	Beige	360-6555
76658-09	Large	Short	Beige	360-6464
76658-03	Large	Regular	Beige	360-6514
76658-15	Large	Tall	Beige	360-6563
76658-10	Extra Large	Short	Beige	360-6472
76658-04	Extra Large	Regular	Beige	360-6522
76658-16	Extra Large	Tall	Beige	360-6571
<b>Legpiece</b>				
76657-00	Extra Small	Regular	Beige	360-6589
76657-05	Extra Small	Tall	Beige	360-6639
76657-01	Small	Regular	Beige	360-6597
76657-06	Small	Tall	Beige	360-6647
76657-02	Medium	Regular	Beige	360-6605
76657-07	Medium	Tall	Beige	360-6654
76657-03	Large	Regular	Beige	360-6613
76657-08	Large	Tall	Beige	360-6662
76657-04	Extra Large	Regular	Beige	360-6621
76657-09	Extra Large	Tall	Beige	360-6670
<b>Footpiece</b>				
76656-05	Extra Small	Regular	Beige	360-6688
76656-00	Extra Small	Long	Beige	360-6738
76656-06	Small	Regular	Beige	360-6696
76656-01	Small	Long	Beige	360-6746
76656-07	Medium	Regular	Beige	360-6704
76656-02	Medium	Long	Beige	360-6753
76656-08	Large	Regular	Beige	360-6712
76656-03	Large	Long	Beige	360-6761
76656-09	Extra Large	Regular	Beige	360-6720
76656-04	Extra Large	Long	Beige	360-6779

## JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732



Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance

### Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for **codefinder**, **jobst.com** on your mobile device.



SCAN ME

**HC21**  
Healthcare

info@hc21.group  
00 353 214 875 055  
healthcare21.eu