

JOBST® FarrowWrap® A-D Below Knee

Ready-to-wear order form and prescription request

For quotes please email: quotes@hc21.group


To order please email Customer Services on: info@hc21.group

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____ Invoice Address: _____

Legpiece



Measuring instructions

Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **A** (floor) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-XL**) & length (**Regular or Tall**) using chart below

A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap legpiece

Footpiece

1. Measure **A1** circumference
2. Measure length **X** from heel to base of little toe
3. Determine size (**XS-XL**) & length (**Regular or Long**) using chart below

JOBST FarrowHybrid (Liners)


1. JOBST FarrowHybrid (Liners) can be used instead of a footpiece, providing calf circumference does not exceed 60cm
2. Measure **A1** circumference
3. Record **shoe size**
4. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



Footpiece



Left:	Shoe Size
Right:	Shoe Size

Measurement only required when ordering JOBST FarrowHybrid (Liners)

	X Small	Small	Medium	Large	X Large
Legpiece					
C Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
B Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
A-D Regular	33-37cm	35-39cm	37-41cm	39-43cm	39-43cm
A-D Tall	38-41cm	40-43cm	42-45cm	44-47cm	44-47cm
NOTE: for A-D length measure posteriorly (back of leg) and follow contours					
Footpiece					
A1 Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
X Regular	16-17cm	17-18cm	19-20cm	20-21cm	22-23cm
X Long	18-19cm	19-20cm	21-22cm	22-23cm	24-25cm
NOTE: for X length measure from the heel to the base of the little toe					

Note: 3cm allowance on either side of measurement

Legpiece

Size: X Small Small Medium Large X Large

Length **A-D** Regular Tall

Fabric: Lite Strong Classic

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

Footpiece

Size: X Small Small Medium Large X Large

Length **X** Regular Long

Fabric: Lite Strong Classic

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST® FarrowHybrid® (Liners)

Size: Small (2-5) Medium (6-8) Large (9-11) X Large (12-13) (shoe size)

Width **A1** Regular (20-27cm) Wide (28-35cm)

Colour: Tan Black (available in Wide variant only)

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST FarrowWrap Order Codes

JOBST® FarrowWrap® Lite (20-30mmHg)

For patients with **mild to moderate** oedema caused by venous or lymphatic conditions, where reduced compression is required

Product Code	Size	Length	Colour	PIP Code
Legpiece				
76654-00	Extra Small	Regular	Tan	360-6233
76654-05	Extra Small	Tall	Tan	360-6282
76654-01	Small	Regular	Tan	360-6241
76654-06	Small	Tall	Tan	360-6290
76654-02	Medium	Regular	Tan	360-6258
76654-07	Medium	Tall	Tan	360-6308
76654-03	Large	Regular	Tan	360-6266
76654-08	Large	Tall	Tan	360-6316
76654-04	Extra Large	Regular	Tan	360-6274
76654-09	Extra Large	Tall	Tan	360-6324
Footpiece				
76652-05	Extra Small	Regular	Tan	360-6332
76652-00	Extra Small	Long	Tan	360-6381
76652-06	Small	Regular	Tan	360-6340
76652-01	Small	Long	Tan	360-6399
76652-07	Medium	Regular	Tan	360-6357
76652-02	Medium	Long	Tan	360-6407
76652-08	Large	Regular	Tan	360-6365
76652-03	Large	Long	Tan	360-6415
76652-09	Extra Large	Regular	Tan	360-6373
76652-04	Extra Large	Long	Tan	360-6423



JOBST® FarrowWrap® Strong (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for skin folds and shape distortion

Product Code	Size	Length	Colour	PIP Code
Legpiece				
76648-00	Extra Small	Regular	Tan	382-0073
76648-05	Extra Small	Tall	Tan	382-0024
76648-01	Small	Regular	Tan	382-0107
76648-06	Small	Tall	Tan	382-0057
76648-02	Medium	Regular	Tan	382-0099
76648-07	Medium	Tall	Tan	382-0040
76648-03	Large	Regular	Tan	382-0081
76648-08	Large	Tall	Tan	382-0032
76648-04	Extra Large	Regular	Tan	382-0065
76648-09	Extra Large	Tall	Tan	382-0016
Footpiece				
76646-04	Extra Small	Regular	Tan	381-9976
76646-09	Extra Small	Long	Tan	381-9927
76646-05	Small	Regular	Tan	382-0008
76646-00	Small	Long	Tan	381-9950
76646-06	Medium	Regular	Tan	381-9992
76646-01	Medium	Long	Tan	381-9943
76646-07	Large	Regular	Tan	381-9984
76646-02	Large	Long	Tan	381-9935
76646-08	Extra Large	Regular	Tan	381-9968
76646-03	Extra Large	Long	Tan	381-9919



JOBST® FarrowWrap® Classic (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for irregular shaped limbs and deep skin folds

Product Code	Size	Length	Colour	PIP Code
Legpiece				
76657-00	Extra Small	Regular	Tan	360-6589
76657-05	Extra Small	Tall	Tan	360-6639
76657-01	Small	Regular	Tan	360-6597
76657-06	Small	Tall	Tan	360-6647
76657-02	Medium	Regular	Tan	360-6605
76657-07	Medium	Tall	Tan	360-6654
76657-03	Large	Regular	Tan	360-6613
76657-08	Large	Tall	Tan	360-6662
76657-04	Extra Large	Regular	Tan	360-6621
76657-09	Extra Large	Tall	Tan	360-6670
Footpiece				
76656-05	Extra Small	Regular	Tan	360-6688
76656-00	Extra Small	Long	Tan	360-6738
76656-06	Small	Regular	Tan	360-6696
76656-01	Small	Long	Tan	360-6746
76656-07	Medium	Regular	Tan	360-6704
76656-02	Medium	Long	Tan	360-6753
76656-08	Large	Regular	Tan	360-6712
76656-03	Large	Long	Tan	360-6761
76656-09	Extra Large	Regular	Tan	360-6720
76656-04	Extra Large	Long	Tan	360-6779

JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732



Do not apply JOBST FarrowHybrid (Liners) if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance

Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for **codefinder**. jobst.com on your mobile device.



info@hc21.group
00 353 214 875 055
healthcare21.eu



JOBST®,
an Essity brand



Essity
T/A BSN medical Limited
PO Box 258 • Willerby • Hull • HU10 6WT

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018 and the (EU) General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit www.healthcare21.eu