

JOBST® Ready-to-Wear

Ready-to-wear order form and prescription request

For quotes please email: quotes@hc21.group

To order please email Customer Services on: info@hc21.group

Date: _____ Purchase Order No.: _____

Measured By: _____

Delivery Address: _____

Patient Name: _____

Tel: _____ Email: _____

Invoice Address: _____

1. <input type="checkbox"/> JOBST® Opaque	<input type="checkbox"/> JOBST® UltraSheer	<input type="checkbox"/> JOBST® Classic	<input type="checkbox"/> JOBST® ForMen Ambition <input type="checkbox"/> JOBST® ForMen Explore
2. <input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights	<input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights	<input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights	<input type="checkbox"/> Knee High
3. RAL Compression Class <input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)	RAL Compression Class <input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)	RAL Compression Class <input type="checkbox"/> CCL 2 (23-32mmHg) <input type="checkbox"/> CCL 3 (34-46mmHg)	RAL Compression Class <input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)
4. <input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe	<input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe	<input type="checkbox"/> Open Toe	<input type="checkbox"/> Closed Toe
5. Silicone Top Band (Thigh High only) <input type="checkbox"/> Sensitive <input type="checkbox"/> Sensitive Wide <input type="checkbox"/> Lace <input type="checkbox"/> Lace Wide <input type="checkbox"/> Dotted <input type="checkbox"/> Dotted Wide	Silicone Top Band (Thigh High only) <input type="checkbox"/> Sensitive <input type="checkbox"/> Sensitive Wide <input type="checkbox"/> Lace <input type="checkbox"/> Lace Wide <input type="checkbox"/> Dotted <input type="checkbox"/> Dotted Wide	Silicone Top Band (Thigh High only) <input type="checkbox"/> Dotted	
6. <input type="checkbox"/> Petite <input type="checkbox"/> Regular	<input type="checkbox"/> Petite <input type="checkbox"/> Regular	<input type="checkbox"/> Petite <input type="checkbox"/> Regular	<input type="checkbox"/> Regular <input type="checkbox"/> Long
7. <input type="checkbox"/> Black <input type="checkbox"/> Natural <input type="checkbox"/> Caramel <input type="checkbox"/> Bronze <input type="checkbox"/> Navy	<input type="checkbox"/> Black <input type="checkbox"/> Natural <input type="checkbox"/> Caramel <input type="checkbox"/> Bronze <input type="checkbox"/> Espresso <input type="checkbox"/> Anthracite	<input type="checkbox"/> Pearl <input type="checkbox"/> Sand	<input type="checkbox"/> Black <input type="checkbox"/> Dark Grey <input type="checkbox"/> Navy <input type="checkbox"/> Khaki <input type="checkbox"/> Brown (Ambition only)
8. <input type="checkbox"/> Size I <input type="checkbox"/> Size II <input type="checkbox"/> Size III <input type="checkbox"/> Size IV <input type="checkbox"/> Size V <input type="checkbox"/> Size VI			

Product Code: -


Quantity Required:

Repeat prescription required every months

Find ordering and prescribing information fast with JOBST Code Finder


No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for codefinder.jobst.com on your mobile device.



Measuring Guide


JOBST® Opaque / JOBST® UltraSheer / JOBST® Classic*



Measurement point	Sizes						Measurements	
Circumference (cm)	I	II	III	IV	V	VI	Left leg	Right leg
cG (no silicone band)	43 - 57cm	45 - 62cm	49 - 67cm	53 - 72cm	56 - 77cm	60 - 81cm		
cG (silicone band wide)	49 - 57cm	53 - 62cm	57 - 67cm	61 - 72cm	65 - 77cm	69 - 81cm		
cG (silicone band)	43 - 48cm	45 - 52cm	49 - 56cm	53 - 60cm	56 - 64cm	60 - 68cm		
cF	39 - 52cm	41 - 56cm	44 - 60cm	47 - 65cm	50 - 69cm	53 - 73cm		
cE	30 - 37cm	33 - 40cm	35 - 43cm	37 - 45cm	39 - 48cm	41 - 51cm		
cD	27 - 33cm	29 - 36cm	32 - 39cm	34 - 42cm	36 - 45cm	38 - 48cm		
cC	28 - 34cm	30 - 37cm	33 - 40cm	35 - 43cm	37 - 46cm	39 - 49cm		
cB1	23 - 27cm	24 - 29cm	26 - 32cm	29 - 35cm	31 - 37cm	33 - 39cm		
cB	18 - 20cm	20 - 22cm	22 - 24cm	24 - 26cm	26 - 28cm	28 - 30cm		
cY	26 - 31cm	28 - 33cm	29 - 35cm	31 - 37cm	32 - 38cm	33 - 40cm		
cA	17 - 22cm	19 - 24cm	21 - 26cm	23 - 29cm	25 - 32cm	27 - 34cm		
Length (cm)	Petite		Regular			Left leg	Right leg	
IK (Tights)	65 - 75cm		75 - 85cm					
IG (Thigh high)	60 - 69cm		70 - 80cm					
ID (Knee high)	34 - 39cm		39 - 45cm					

*JOBST Classic thigh high garment available in Regular length only

JOBST® ForMen Ambition / JOBST® ForMen Explore



Measurement point	Sizes						Measurements	
Circumference (cm)	I	II	III	IV	V	VI	Left leg	Right leg
cD	28 - 34cm	31 - 38cm	34 - 41cm	37 - 46cm	40 - 49cm	42 - 52cm		
cC	30 - 36cm	32 - 40cm	34 - 43cm	38 - 47cm	40 - 50cm	43 - 53cm		
cB1	23.5 - 27.5cm	26 - 30.5cm	28.5 - 33cm	31 - 36.5cm	33.5 - 39cm	35.5 - 41.5cm		
cB	18 - 20cm	20 - 22cm	22 - 24cm	24 - 26cm	26 - 28cm	28 - 30cm		
cY	26 - 31cm	28 - 33cm	29 - 35cm	31 - 37cm	32 - 38cm	33 - 40cm		
cA	17 - 22cm	19 - 24cm	21 - 26cm	23 - 29cm	25 - 32cm	27 - 34cm		
Length (cm)	Regular			Long			Left leg	Right leg
ID (Knee high)	39 - 45cm			45 - 51cm				

Please note

Some patients may fit into more than one size, use clinical judgement to select most appropriate size. If the best fit is not clear, then the patient may require a made to measure flat knit garment or an adjustable compression wrap device.

Measuring instructions

1. Measure each circumference point for the garment style required (e.g. knee high, thigh high or tights). Record each measurement in the left & right leg columns
2. Starting with cA, circle each measurement range along the row that the circumference measurement fits into. Repeat for each circumference point needed for the required garment style (see example 1 below)
3. Determine the size required based on the size column with the most circles (see example 2 below)
4. Measure the leg length, taking a straight measurement from the floor to the required measuring point (do not follow the contour of the leg). Record length measurements in the left & right leg columns. Determine the appropriate leg length (Petite or Regular) or (Regular or Long)

1 Example: If cA measurement is 22cm

Circumference (cm)	I	II	III
cG (tights only)	43 - 57cm	45 - 62cm	49 - 67cm
cG (silicone band wide)	49 - 57cm	53 - 62cm	57 - 67cm
cG (silicone band)	43 - 48cm	45 - 52cm	49 - 56cm
cF	39 - 52cm	41 - 56cm	44 - 60cm
cE	30 - 37cm	33 - 40cm	35 - 43cm
cD	27 - 33cm	29 - 36cm	32 - 39cm
cC	28 - 34cm	30 - 37cm	33 - 40cm
cB1	23 - 27cm	24 - 29cm	26 - 32cm
cB	18 - 20cm	20 - 22cm	22 - 24cm
cY	26 - 31cm	28 - 33cm	29 - 35cm
cA	17 - 22cm	19 - 24cm	21 - 26cm

2 Size III is the most appropriate size

	Sizes						Meas
	II	III	IV	V	VI	Left	
cG (tights only)	45 - 62cm	49 - 67cm	53 - 72cm	56 - 77cm	60 - 81cm		
cG (silicone band wide)	53 - 62cm	57 - 67cm	61 - 72cm	65 - 77cm	69 - 81cm		
cG (silicone band)	45 - 52cm	49 - 56cm	53 - 60cm	56 - 64cm	60 - 68cm		
cF	41 - 56cm	44 - 60cm	47 - 65cm	50 - 69cm	53 - 73cm		
cE	33 - 40cm	35 - 43cm	37 - 45cm	39 - 48cm	41 - 51cm		
cD	30 - 36cm	32 - 39cm	34 - 42cm	36 - 45cm	38 - 48cm		34
cC	30 - 37cm	33 - 40cm	35 - 43cm	37 - 46cm	39 - 49cm		35
cB1	24 - 29cm	26 - 32cm	29 - 35cm	31 - 37cm	33 - 39cm		29
cB	20 - 22cm	22 - 24cm	24 - 26cm	26 - 28cm	28 - 30cm		23
cY	28 - 33cm	29 - 35cm	31 - 37cm	32 - 38cm	33 - 40cm		33
cA	19 - 24cm	21 - 26cm	23 - 29cm	25 - 32cm	27 - 34cm		22

Measuring Points	
Circumference	Length
cG = widest part of thigh	IK = floor to crotch
cF = middle of thigh	IG = floor to gluteal fold
cE = middle of patella with leg flexed at 30-45 degrees	ID = floor to two fingers width below crease of the knee
cD = two fingers width below crease of the knee	
cC = widest part of calf	
cB1 = insertion of the Achilles tendon	
cB = narrowest part of ankle	
cY = heel to front ankle crease (maximum dorsiflexion)	
cA = base of toes	

Check out our measuring guide videos

Scan the relevant QR code for help with measuring for JOBST® Ready-to-Wear (RtW)



Knee High



Thigh High

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, and the EU General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit [Group Data Protection Statement - Healthcare21](#)

E: info@hc21.group

W: www.healthcare21.eu



info@hc21.group
00 353 214 875 055
healthcare21.eu

Customer Service: 0818 21 21 21

Compression Therapy Helpline: 0800 389 8424



JOBST®, an Essity brand

