

JOBST® UlcerCARE

2-in-1 Compression System

Ready-to-wear order form and prescription request

For quotes please email: quotes@hc21.group

To order please email Customer Services on: info@hc21.group

Date: _____ Purchase Order No.: _____


Patient Name: _____

Measured By: _____

Tel: _____ Email: _____

Delivery Address: _____

Invoice Address: _____



Measuring Instructions

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Determine size (**S-4XL**) using chart below

Size	S	M	L	XL	2XL	3XL	4XL
C	29 - 35cm	34 - 41cm	40 - 47cm	46 - 54cm	48 - 56cm	50 - 58cm	52 - 60cm
B	18 - 21cm	21 - 24cm	24 - 27cm	27 - 30cm	30 - 33cm	33 - 36cm	36 - 39cm

Top Tip

If the measurements fall into two sizes (e.g. B is Medium and C is Large), choose the size with the measurement that is closest to the middle of the range.

For example, B is 22cm and C is 42cm.

22cm is in the middle of the Medium size (21-24cm).

42cm is at the start of the Large size (40-47cm).

In this example, we would advise you order the Medium garment.

Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST® UlcerCARE



JOBST® UlcerCARE Kit

Each kit contains one outer stocking and two liners

Size

- Small Medium Large
 X Large 2X Large 3X Large
 4X Large

Colour

- Beige Black

Zipper (to be placed on the side opposite to the ulceration)

- No zip Zip left Zip right

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST® UlcerCARE Liners

Each pack contains three liners

Size

- Small Medium Large
 X Large 2X Large 3X Large
 4X Large

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST® UlcerCARE Order Codes

JOBST® UlcerCARE Kit (40mmHg)

- 1 x Outer Stocking
- 2 x Liners
- With Zip

Product Code	Size
Zipper to be placed on the side opposite to the ulceration	
Zipper Left - Black	
73631-41	S
73631-42	M
73631-43	L
73631-44	XL
73631-45	2XL
73631-46	3XL
73631-47	4XL
Zipper Right - Black	
73631-71	S
73631-72	M
73631-73	L
73631-74	XL
73631-75	2XL
73631-76	3XL
73631-77	4XL
Zipper Left - Beige	
73631-21	S
73631-22	M
73631-23	L
73631-24	XL
73631-25	2XL
73631-26	3XL
73631-27	4XL
Zipper Right - Beige	
73631-61	S
73631-62	M
73631-63	L
73631-64	XL
73631-65	2XL
73631-66	3XL
73631-67	4XL

JOBST® UlcerCARE Kit (40mmHg)

- 1 x Outer Stocking
- 2 x Liners
- Without Zip

Product Code	Size	PCRS Code
Black		
73630-41	S	86342
73630-42	M	86342
73630-43	L	86342
73630-44	XL	86342
73630-45	2XL	86342
73630-46	3XL	86342
73630-47	4XL	86342
Beige		
73630-21	S	86342
73630-22	M	86342
73630-23	L	86342
73630-24	XL	86342
73630-25	2XL	86342
73630-26	3XL	86342
73630-27	4XL	86342

JOBST® UlcerCARE Liners (15-20mmHg)

- 3 x Liners

Product Code	Size	PCRS Code
Replacement Liners - White		
73632-21	S	86084
73632-22	M	86084
73632-23	L	86084
73632-24	XL	86084
73632-25	2XL	86084
73632-26	3XL	86084
73632-27	4XL	86084

Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for codefinder.jobst.com on your mobile device.



Scan to access HC21 Compression Products available on PCRS



	info@hc21.group
	00 353 214 875 055 healthcare21.eu