

By your side
and beyond




Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set

Discover the future of enteral feeding

The **Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set** has been designed to deliver thick formula. Thick formulas are defined as enteral fluids of smooth consistency that would be categorized as level 2, 3 or 4 drinks within the International Dysphagia Diet Standardisation Initiative (IDDSI) framework. These versatile sets are not made with DEHP, latex or BPA.



Explore the advanced capabilities of Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Sets



Hydration

Maintain tube patency and hydration when used with the Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag

Innovative

Delivers thick formulas

Easy to use

Designed for use with ENFit™ connections and color-coded, patented cassettes that load with one hand



Ordering information

Code	Description	Volume	Sales Unit
E5FDR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set, Non-Sterile	500 mL	30
E10FDR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set Non-Sterile	1000 mL	30
E5FFR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag, Non-Sterile	500 mL	30
E10FFR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag, Non-Sterile	1000 mL	30
E5FDSR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set, - Sterile	500 mL	30
E10FDSR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set, - Sterile	1000 mL	30
E10FFSR	Kangaroo OMNI™ ENTelliSet™ Thick Formula Feeding Set with Flush Bag, - Sterile	1000 mL	30

For more information or to order, contact your sales representative or visit cardinalhealth.co.uk

*Confirm with your local Cardinal Health Care representative the availability of Non-Sterile Kangaroo OMNI™ Enteral Feeding Sets in your country.

For healthcare professionals only.

© 2024 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, KANGAROO, the Kangaroo LOGO, KANGAROO OMNI and KANGAROO ENTELLISET are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent cardinalhealth.com/patents. Lit. No. 3MS23-2633760 (05/2024)

