



User height  
850–1700 mm



Max. load  
60 kg



Environment  
Indoor & outdoor

# Kudu®

The Kudu is a practical wheelchair which combines form and function to provide the right fit and support for a comfortable posture.

The Kudu is available in 4 sizes.

## Highlights about Kudu



### Growth adaptable

The base of the seat and the backrest are equipped with a mechanism that makes it possible to adjust height and depth as well as mount fix-locks for vests and belts.



Accessories are easily added and removed as needed for maximum user comfort



### Centre of gravity

To adapt to different needs, the seat has a centre of gravity making it possible to adjust the rear wheels in 3 different positions.



### Safe transportation

The Kudu complies with ISO 7171-19 for transporting children up to 57 kg in motor vehicles.



## Standard

From -3 to 45° adjustable tilt, -5 to -30 back recline, height and angle adjustable foot supports, seat depth growth of 100 mm, poke protectors.

## Options

### Sizes

- Size 1
- Size 2
- Size 3
- Size 4

### Cushions

- Kudu cushions, upcharge
- No cushions
- Cushions with contour, upcharge

### For the hands

- Push brace, angle adjustable
- Push brace, height and angle adjustable
- Push handles, height adjustables

### Foot support

- Foot plate, centrally mounted
- Parted foot rests (available for size 2+3+4)

### Rear wheels

- 12½"/16" solid or air
- 20"/22"/24" solid or air with quick release

### Front wheels

- 6" solid front wheels
- 7" solid or air front wheels

### Brakes

- Wheel brakes (available for 12½", 20", 22" and 24")
- Drum brakes

### Anti-tips & spoke protectors

- Anti-tips
- Spoke protectors

### Frame colour



Anthracite grey



White



Red
















Kudu with Push brace, angle adjustable



Kudu with Push handles, height adjustable

## Measurements

		1	2	3	4
	Effective seat width	260 mm	300 mm	340 mm	380 mm
	Effective seat depth	220–320 mm	270–370 mm	310–410 mm	350–450 mm
	Backrest height	370–460 mm	400–520 mm	430–580 mm	430–630 mm
	Seat surface height at front edge	480 mm	480 mm	480 mm	480 mm
	Max. load	40 kg	50 kg	60 kg	60 kg
	Max. load at transportation	40 kg	50 kg	57 kg	57 kg
	Backrest angle	-5° - 30°	-5° - 30°	-5° - 30°	-5° - 30°
	Seat plane angle	-3° - 45°	-3° - 45°	-3° - 45°	-3° - 45°
	Length	690 mm	690 mm	740 mm	740 mm
	Overall width	590 mm	630 mm	670 mm	710 mm
	Folded length	825 mm	825 mm	875 mm	875 mm
	Folded height	625 mm	625 mm	625 mm	625 mm
	Total mass	19 kg	21 kg	23 kg	25 kg

The Kudu is measured with the following wheel combination, placed in the standard wheel position: 22" rear wheel and 6" front wheel. Further measurements at our website

## Nice to know

The many different configuration possibilities and accessories for the Kudu makes it ideal for situations where complex and individualized seating solutions is needed, also when reusing of the product is desirable. Examples of individualized solutions are, placing the centrally mounted foot plate underneath the seat or using the angle adjustable armrest and tray, allowing the tray to be horizontal when the seat plate is angled.

The focus of the seating solution will often be comfort, however when configured with large driving wheel and cambering it is possible to promote self-propelling of the wheelchair over short distances to support user independence.



## Accessories



**Head support**

99xxxx-xx

\*see accessories category



**Back extension cushion**

95743-xxx



**Fixed side supports**

95401-xx



**Swing-away side supports**

8310xx-xxx

\*see accessories category



**Back wedges**

E6301x-60



**Arm rests/Tray**

95712-x/9576351



**Hip supports**

95713-x0



**Pommel**

95430-x0

**Padding**

95712-x

**Padding**

95430-x0260



**Vest**

89113-x



**Mud guards**

9520xx



**Transport fittings**

95900-16



**Adductors**

95429-10

For more accessories - visit our website





**Vent tray**  
95281-10



**Pommel**  
95430-x



**Upper arm supports**  
8910407-x0



**Angle adjustable arm rests, incl. side plates**  
95712-x0



**Canopy**  
87401x0  
(Requires 9589100)



**Thermo bag**  
87460-xxx



**Rain cape with hood**  
87465-xxx

