

Levamed

Ankle support with silicone cushions




Indications

- slight instability of the ankle joint
- sprains and distorsions
- joint effusions and swellings due to arthritis and osteoarthritis
- chronic, postoperative and posttraumatic inflammation

Mode of action

- massage effect due to integrated silicone pads
→ reduction of swellings (oedema, haematomas)
- proprioceptive compression
→ guidance and stabilisation of the joints

Product benefits



Innovative comfort knit with anatomically tailored support, soft instep and supple edging.




Integrated silicone pad with a 3D profile for a gentle massage effect and to reduce swelling.

Accessories



A donning aid / slip-in sock is supplied with each product.



silver




sand



Clima Comfort



Clima Fresh



Comfort Zone



medi Kidz – in sizes I + II also for children suitable.



circumference ankle (cm)	17 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30
size	I	II	III	IV	V	VI
article number	K.010.Δ11	K.010.Δ12	K.010.Δ13	K.010.Δ14	K.010.Δ15	K.010.Δ16

Δ = fill in for colour: 2 = silver, 8 = sand