

# iVascular<sup>®</sup>

therapies for living


## luminor 14m

Paclitaxel eluting PTA balloon dilatation catheter

**Safe and precise**  
delivery of Paclitaxel

**LUMINOR Registry**  
207 patients

Nanotechnology  
dosage system.  
**Ultra-thin** coating



Assured coating  
**stability** during  
navigation



Proprietary coating  
technology **TransferTech**



Homogeneous  
& **anti-cracking** coating



[www.ivascular.global](http://www.ivascular.global)

# luminor 14m features

- > **Over-the-wire catheter (OTW)**
- > Usable catheter length: **100 cm** and **150 cm**
- > Recommended guidewire: **0.014"**
- > Tip profile: **0.017"**
- > Crossing profile: **from 0.026" up to 0.042"**
- > Nominal pressure: **7 atm**
- > RBP: **16 atm**
- > ABP: **22 atm**

Product with CE mark, certified by Notified Body 0318

Usable catheter length	Balloon Ø (mm)	Balloon Length (mm)					
		40	60	80	120	150	200
100 cm	1.5	BP DPC14M 100 150 040	BP DPC14M 100 150 060	BP DPC14M 100 150 080	-	-	-
	2.0	BP DPC14M 100 200 040	BP DPC14M 100 200 060	BP DPC14M 100 200 080	BP DPC14M 100 200 120	BP DPC14M 100 200 150	BP DPC14M 100 200 200
	2.5	BP DPC14M 100 250 040	BP DPC14M 100 250 060	BP DPC14M 100 250 080	BP DPC14M 100 250 120	BP DPC14M 100 250 150	BP DPC14M 100 250 200
	3.0	BP DPC14M 100 300 040	BP DPC14M 100 300 060	BP DPC14M 100 300 080	BP DPC14M 100 300 120	BP DPC14M 100 300 150	BP DPC14M 100 300 200
	3.5	BP DPC14M 100 350 040	BP DPC14M 100 350 060	BP DPC14M 100 350 080	BP DPC14M 100 350 120	BP DPC14M 100 350 150	BP DPC14M 100 350 200
	4.0	BP DPC14M 100 400 040	BP DPC14M 100 400 060	BP DPC14M 100 400 080	BP DPC14M 100 400 120	-	-
150 cm	1.5	BP DPC14M 150 150 040	BP DPC14M 150 150 060	BP DPC14M 150 150 080	-	-	-
	2.0	BP DPC14M 150 200 040	BP DPC14M 150 200 060	BP DPC14M 150 200 080	BP DPC14M 150 200 120	BP DPC14M 150 200 150	BP DPC14M 150 200 200
	2.5	BP DPC14M 150 250 040	BP DPC14M 150 250 060	BP DPC14M 150 250 080	BP DPC14M 150 250 120	BP DPC14M 150 250 150	BP DPC14M 150 250 200
	3.0	BP DPC14M 150 300 040	BP DPC14M 150 300 060	BP DPC14M 150 300 080	BP DPC14M 150 300 120	BP DPC14M 150 300 150	BP DPC14M 150 300 200
	3.5	BP DPC14M 150 350 040	BP DPC14M 150 350 060	BP DPC14M 150 350 080	BP DPC14M 150 350 120	BP DPC14M 150 350 150	BP DPC14M 150 350 200
	4.0	BP DPC14M 150 400 040	BP DPC14M 150 400 060	BP DPC14M 150 400 080	BP DPC14M 150 400 120	-	-

● Introducer compatibility: **4F**

Distributed by:




**Aquilant**

Aquilant House, Unit B1-B2 Bond Close,  
Kingsland Business Park, Basingstoke,  
Hampshire, RG248PZ, United Kingdom

T: +44 (0) 1256 365 454  
F: +44 (0) 1256 365 486

E: [contactus@aquilantservices.com](mailto:contactus@aquilantservices.com)  
W: [www.aquilantservices.com](http://www.aquilantservices.com)



Manufactured by:

**Life Vascular Devices Biotech, S.L.**  
[www.ivascular.global](http://www.ivascular.global)  
[info@ivascular.global](mailto:info@ivascular.global)