

# iVascular<sup>®</sup>

therapies for living

## luminor 35


Paclitaxel eluting PTA balloon dilatation catheter

**EFFPAC trial**  
171 patients



Homogeneous  
& **anti-cracking**  
coating

**Safe and precise**  
delivery of Paclitaxel



Assured coating  
**stability** during  
navigation

Nanotechnology  
dosage system.  
**Ultra-thin** coating



Proprietary coating  
technology **TransferTech**



[www.ivascular.global](http://www.ivascular.global)

# luminor 35 features

- > **Over-the-wire catheter (OTW)**
- > Usable catheter length: **80 cm** or **140 cm**
- > Recommended guidewire: **0.035"**
- > Tip profile: **0.036"**
- > Crossing profile: **from 0.065" up to 0.077"**
- > Nominal pressure: **6/7 atm**
- > RBP: **16 atm**
- > ABP: **23 atm**

Product with CE mark, certified by Notified Body 0318

Usable catheter length	Balloon Ø (mm)	Balloon Length (mm)					
		20	40	60	80	120	150
80 cm	5.0	BP DPC35 080 500 020	BP DPC35 080 500 040	BP DPC35 080 500 060	BP DPC35 080 500 080	BP DPC35 080 500 120	BP DPC35 080 500 150
	6.0	BP DPC35 080 600 020	BP DPC35 080 600 040	BP DPC35 080 600 060	BP DPC35 080 600 080	BP DPC35 080 600 120	BP DPC35 080 600 150
	7.0	BP DPC35 080 700 020	BP DPC35 080 700 040	BP DPC35 080 700 060	BP DPC35 080 700 080	BP DPC35 080 700 120	-
140 cm	5.0	BP DPC35 140 500 020	BP DPC35 140 500 040	BP DPC35 140 500 060	BP DPC35 140 500 080	BP DPC35 140 500 120	BP DPC35 140 500 150
	6.0	BP DPC35 140 600 020	BP DPC35 140 600 040	BP DPC35 140 600 060	BP DPC35 140 600 080	BP DPC35 140 600 120	BP DPC35 140 600 150
	7.0	BP DPC35 140 700 020	BP DPC35 140 700 040	BP DPC35 140 700 060	BP DPC35 140 700 080	BP DPC35 140 700 120	-

- Introducer compatibility: **5F**
- Introducer compatibility: **6F (ID 2.20 mm)**

Distributed by:



**Aquilant**

Aquilant House, Unit B1-B2 Bond Close,  
Kingsland Business Park, Basingstoke,  
Hampshire, RG248PZ, United Kingdom

T: +44 (0) 1256 365 454  
F: +44 (0) 1256 365 486

E: [contactus@aquilantservices.com](mailto:contactus@aquilantservices.com)  
W: [www.aquilantservices.com](http://www.aquilantservices.com)



Manufactured by:

Life Vascular Devices Biotech, S.L.  
[www.ivascular.global](http://www.ivascular.global)  
[info@ivascular.global](mailto:info@ivascular.global)