



"WITHIN 1 HOUR OF DELIVERY I WAS IN THE BUSH AND ON THE BEACH! A DREAM COME TRUE AND PROBABLY ONE OF THE BEST FEELINGS I'VE EVER HAD, ABSOLUTE FREEDOM."

Nick
New South Wales

magic
Mobility™

Extreme X8

**OFF-ROAD
AND OFF
THE MAP**

HIT THE ROAD

Or the sand, or the bush, or the snow.*

*Not to be driven through water

MAKING YOUR CHAIR WORK FOR YOU

Our wheelchairs are found on building sites and farms, in stables and national parks - allowing people to continue to work.



DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help challenge you to embrace your freedom.



V7/AUST/FEB22

TAILORED FOR
YOUR LIFESTYLE



AUSTRALIAN
MADE

Extreme
X8™



MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. Sand, snow and every terrain you can throw at it. With unrivalled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

"WE'VE BEEN
CREATING
MADE-TO-ORDER
WHEELCHAIRS
FOR OVER 25
YEARS, SO THAT
YOU CAN LIVE THE
LIFE YOU
CHOOSE"

Ashley Daff
R&D Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

Basic specifications:*

- Length: 1030mm
- Width: 710mm
- Seat to floor height: 470mm (490mm with swing away leg rests)
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Gel batteries (70Ah)
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide).

Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power or gas elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

Options and accessories:

- Battery upgrade 90Ah (provides 20% more battery range)
- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Puncture resistant tyre liners
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- Crash test compliant tie downs
- RAM® mounting systems for mobile phones, tablets, cameras and cup holders.

The right wheelchair for:**

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living.

Base options:



FLAME RED - BLACK RIMS



OCEAN BLUE - SILVER RIMS



ARMY CAMO - SILVER RIMS

*Specifications subject to change. Check online for current version
**Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au
Explore with us magicmobility.com.au



Scan for agent locations