

medi 4C[®] flex

Hyperextension brace with button closer or hook-closer

Indications

- stable vertebral compression fractures of the lower thoracic and lumbar spine
- medium-term immobilisation after intervertebral disc surgery
- axial stabilisation in osteoporotic bone collapse of the thoracic and lumbar spine
- temporary orthotic care after open surgery of unstable vertebral body fractures
- permanent orthotic care for inoperable spinal tumors and metastases
- post-operative orthotic care after surgery for tumors and metastases of vertebral bodies
- aftercare of spinal decompression surgery with and without internal fixation

Mode of action

- hyperextension of the vertebral column provides stabilisation in the saggital plane
- rotational stability is restricted by the subclavicular extension pads
- stable framework construction restricts lateral flexion

Features

- moveable hip bar for better comfort while sitting
- available in 4 sizes with optional adjustments of the side bar length
- easy to use with new developed button closer
- width of chest bars is adjustable
- also available as short version and ultrashort version
- available with button or hook closer



medi 4C flex



button closer