

medi Humeral fracture brace

Indications

- Postoperative care of subcapital or shaft fractures to the humerus
- Conservative care of humerus fractures

Contraindications

- None currently known

Product benefits



Rigid bars removable for custom fit and x-ray



Removable upper arm belts

Option to reposition or remove upper arm belt for flexible range of application and compression



Lower arm sleeve


Easy to open for passive exercises

Breathable material ensures maximum comfort



medi Humeral fracture brace

	size	circumference	article number
Right	S	21 - 28 cm	R090312
	M	26 - 34 cm	R090313
	L	31 - 40 cm	R090314
Left	S	21 - 28 cm	R090212
	M	26 - 34 cm	R090213
	L	31 - 40 cm	R090214



measurement:
circumference of biceps



4 104611415067681 97E15 / 03.2013