

medi ROM[®]

Universal knee splint with flexion/extension limitation and options for immobilisation

Indications

- post-operative immobilisation
- for immobilisation after knee injuries and for later mobilisation
- for secure stabilisation of injured or operated collateral and/or cruciate ligaments

Mode of action


- collateral stabilisation of the knee joint
- post-operative immobilisation
- controlled mobilisation

Features

- new, patented “QuickSet” hinge (very easy to use)
- limitation in steps of 10° to 15°
flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
extension: 0°, 10°, 20°, 30°
immobilisation: 0°, 10°, 20°, 30°
- easy, comfortable ‘click closing’ strap fasteners
- malleable, metal side bars
- adjustable skin friendly padding
- also available as “light” version (= medi ROM cool)
- medi ROM / medi ROM cool available as short and long version
- numerous options to shorten the side bars



medi ROM



medi ROM cool



new, patented “QuickSet” hinge

medi ROM / medi ROM cool			
size	length cm	article number	article number
		medi ROM	medi ROM cool
short	57	G.180.009	G.184.111
long	63	G.180.001	G.180.111